

MX-PRO®

bariatric transport

stryker®



EMS Equipment

3800 E. Centre Ave
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 784 4336

www.ems.stryker.com

Standard Features

- Transport at load height capability
- Positive action height adjustment
- Easy to use release handle design
- Color-coded controls
- Dampened action during hot drops
- Lap belt extensions
- High visibility powder-coated frame
- Light-weight, rugged aluminum construction
- Scientifically-optimized lift bar design
- Lower lifting bar
- Seven height positions
- Integrated bumper system
- Lift-capable safety bar
- Perforated litter surface
- One-hand release breakaway head section with safety bar
- Floor-mounted safety hook
- One-hand release fold-down siderails
- One-hand release, infinite positioning, lift assist pneumatic backrest
- Oversized wheels with sealed caster and wheel bearings
- Reflective labeling
- Sealed flat mattress
- Shock position
- Two lap belts and one four-point shoulder restraint

Optional Features

- Heavy duty two or three-stage IV poles (patient right or left)
- Base storage tray
- Non-pocketed head-end storage
- Height limit kit*
- Permanent or removable oxygen bottle holders (foot end)
- Rigid head and foot end push/pull handle
- Single or dual wheel lock
- X-frame guards
- Patient security straps
- Side lift handles
- Tow package
- Transfer-Flat (9' x 5', 1,600 lb capacity)

Specifications

• Model Number	6083	
• Height¹		
Position 1	13.5"	(34 cm)
Position 2	24.5"	(62 cm)
Position 3	27.5"	(70 cm)
Position 4	30"	(76 cm)
Position 5	32.5"	(83 cm)
Position 6	34.5"	(88 cm)
Position 7	37"	(94 cm)
• Length		
Standard	80"	(203 cm)
Minimum	62"	(157 cm)
• Width	29"	(74 cm)
• Weight²	111 lbs	(50 kg)
• Wheels		
Diameter	6"	(15 cm)
Width	2"	(5 cm)
• Articulation		
Backrest	0° - 73°	
Shock Position	+14°	
• Maximum Weight Capacity	1600/850 lbs	(727/386 kg)
• Minimum Operators Required		
Occupied Cot	2	
Unoccupied Cot	1	
• Recommended Fastener System		
Floor Mount	Model 6370 or 6377	
Wall Mount	Model 6371	
• Recommended Loading Height*	Up To 33"	(84 cm)

Warranty

- One year parts, labor and travel or two years parts

*Height limit kit recommended for load heights less than 28" (71 cm).
High height kit recommended for load heights more than 32" (81 cm).

¹Height measured from bottom of mattress at seat section to ground level.

²Cot is weighed without mattress and restraints.

Stryker reserves the right to change specifications without notice.

In-service video included with every order.

The MX-PRO® Bariatric Transport is designed to be compatible with competitive cot fastener systems.