

Power-PRO™

powered ambulance cot

stryker®



Standard Features

- (2) DEWALT® 24 volt industrial grade batteries
- DEWALT® industrial grade battery charger (fan cooled)
- Automatic in-cot fastener shut-off
- Manual back-up system
- Automatic high-speed retract
- Battery powered hydraulic lift system
- Preventative maintenance hour meter
- Battery power indicator
- Settable load height with jog function
- Color-coded controls
- High visibility powder-coated frame
- Light-weight, rugged aluminum construction
- Scientifically-optimized lift bar and control design
- Lower lifting bar
- Lift-capable safety bar
- Integrated bumper system
- Perforated litter surface
- One-hand release retracting head section with safety bar
- Floor-mounted safety hook
- One-hand release, fold-down side rails
- One-hand release, infinite positioning, lift assist pneumatic backrest
- Oversized wheels with sealed caster and wheel bearings
- Sealed bolster mattress
- Shock position
- Two lap belts and one four-point shoulder restraint
- Single wheel lock
- Pull handle
- X-frame guards
- Power washable

Optional Features

- Heavy duty two- or three-stage IV poles (patient right or left)
- Permanent or removable oxygen bottle holders (foot end)
- Sealed flat mattress
- Premium mattress
- Dual wheel locks
- Head extension
- Pillow
- Equipment hook
- Backrest storage pouch
- Head end storage flat
- DEWALT® battery charger (12v, 110v and 220v models available)

EMS Equipment

3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 784 4336

www.ems.stryker.com

Specifications

Model Number	6500	
Height¹ (infinite height positioning between lowest and highest position)		
Lowest Position	14"	(36 cm)
Highest Position	41.5"	(105 cm)
Length		
Standard	80"	(203 cm)
Minimum	63"	(160 cm)
Width	23"	(58 cm)
Weight²	125 lbs	(57 kg)
Wheels		
Diameter	6"	(15 cm)
Width	2"	(5 cm)
Articulation		
Backrest	0 – 73°	
Shock Position	+15°	
Maximum Weight Capacity³	700 lbs	(318 kg)
Minimum Operator Required		
Occupied Cot	2	
Unoccupied Cot	1	
Recommended Fastener System		
Floor Mount	Model 6370 or 6377	
Wall Mount	Model 6371	
Recommended Loading Height⁴	Up to 36"	(91 cm)

Warranty

- One-year parts, labor and travel warranty on all soft goods
- Two-year parts, labor and travel warranty on entire cot excluding soft goods
- Three-year parts, labor and travel warranty on X-frame components
- Three-year limited warranty on powertrain components

¹Height measured from bottom of mattress, at seat section, to ground level.

²Cot is weighed with one battery, without mattress and restraints.

³700 lb weight capacity with an unassisted lift capacity of 500 lbs.

⁴Load wheel height can be set between 26" and 36".

Stryker reserves the right to change specifications without notice.

In-service video included with every order.

The Power-PRO is designed to conform to the Federal Specification for the Star-of-Life Ambulance KKK-A-1822E.

The Power-PRO is designed to be compatible with competitive cot fastener systems.

Patents pending.

©Stryker Corporation
06-SEM-1408
Mkt Lit-152-20 Jul 2006 Rev A



DeWALT® is a registered trademark of Black & Decker Inc.