

Power-PRO™ IT

powered incubator transport

stryker®



EMS Equipment

3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2315
toll free: 800 784 4336

www.ems.stryker.com

Standard Features

- (2) DEWALT® 24 volt industrial grade batteries
- DEWALT® industrial grade battery charger (fan cooled)
- Manual back-up system
- Automatic high-speed retract
- Battery powered hydraulic lift system
- Preventative maintenance hour meter
- Battery level indicator
- Color-coded controls
- High visibility powder-coated frame
- Light-weight, rugged aluminum construction
- Scientifically optimized lift bar and control system
- Integrated bumper system
- Oversized wheels with sealed caster and wheel bearings
- X-Frame Guards
- Power washable
- Tie down brackets
- Compatible with airsled or equivalent configurations

Optional Features

- Rigid push bars
- Litter extension
- Head end storage flat
- Base storage net
- DEWALT® battery charger (12v, 110v, and 220v models available)

Compatibility

- Airborne Stackable
- Airborne Side-by-Side
- Dräger Stackable
- Dräger Extended
- No Adapter Option

Specifications

Model Number	6510
Height ¹ (infinite height positioning between lowest and highest position)	
Highest Position	41.5" (105 cm)
Lowest Position	14" (36 cm)
Length	
Standard	63" (160 cm)
Maximum	81" (206 cm)
Width	23" (58 cm)
Weight ²	134 lbs (61 kg)
Wheels	
Diameter	6" (15 cm)
Width	2" (5 cm)
Max. Weight Capacity	700 lbs (318 kg)
Minimum Operators Required	
Occupied Cot	2
Unoccupied Cot	1
Recommended Fastener System	
Floor Mount	Model 6370 or 6377
Wall Mount	Model 6371
Recommended Loading Height ³	Up to 36" (91 cm)

Warranty

- Two-year complete cot (one-year soft goods)
- Three-year X-Frame components
- Three-year limited powertrain

Extended warranties available

¹ Height measured from the top of cot, at the center point, to ground level.

² Cot is weighed with one battery, without incubator.

³ Can accommodate load decks up to 36". Load wheel height can be set between 26" and 36".

Stryker reserves the right to change specifications without notice.

In-service video included with each order.

The Power-PRO IT is recognized by Underwriters Laboratory to IEC 60601.

The Power-PRO IT is designed to conform to Federal Specification for the Star-of-Life Ambulance KKK-1822-E.

The Power-PRO IT is designed to be compatible with competitive cot fastener systems.

Patents pending.



DEWALT® is a registered trademark of Black & Decker Inc.