

# STAIR-PRO<sup>®</sup>

## stair chairs

**stryker<sup>®</sup>**



### EMS Equipment

3800 E. Centre Ave  
Portage, MI 49002 USA  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 784 4336

www.ems.stryker.com

### Standard Features

- High visibility powder-coated frame
- Color-coded controls
- Extendable head and foot end lift handles
- Molded hand grips
- Light-weight, rugged aluminum construction
- Oversized rear wheels with sealed bearings
- Dual wheel locks
- Compact storage size
- Positive action locking mechanism
- Power-washable
- Grease-free maintenance

### Model 6250 Standard Features

- Conforms to BS EN 1865

### Model 6251 Standard Features

- 4" (10 cm) front caster wheels

### Model 6252 Standard Features

- Innovative Stair-TREAD™ system
- Extendable upper control handle
- 4" (10 cm) front caster wheels

### Optional Features

- Locking head end lift handles
- Two-piece molded ABS seat
- One-piece traditional vinyl seat
- Color-coded patient restraints<sup>4</sup> (plastic buckles)
- Non-absorbent patient restraints<sup>4</sup> (metal buckles)
- Polyester patient restraints<sup>4</sup> (metal buckles)
- Patient head support (model 6252 only)
- O2 bottle holder
- Vinyl head support
- Polypropylene head support

*In-service video included with every order.*

### Specifications<sup>1</sup>

| • Models                    | 6250    |          | 6251    |          | 6252     |          |
|-----------------------------|---------|----------|---------|----------|----------|----------|
| • Height                    | 36"     | (91 cm)  | 36"     | (91 cm)  | 37.5"    | (95 cm)  |
| • Width                     | 20.5"   | (52 cm)  | 20.5"   | (52 cm)  | 20.5"    | (52 cm)  |
| • Depth <sup>2</sup>        | 28"     | (71 cm)  | 28"     | (71 cm)  | 28"      | (71 cm)  |
| • Folded Depth              | 8"      | (20 cm)  | 8"      | (20 cm)  | 8"       | (20 cm)  |
| • Weight                    | 20 lbs  | (9 kg)   | 23 lbs  | (10 kg)  | 31.5 lbs | (14 kg)  |
| • Maximum Load <sup>3</sup> | 500 lbs | (227 kg) | 500 lbs | (227 kg) | 500 lbs  | (227 kg) |

### Warranty

- One year parts, labor and travel or two years parts

<sup>1</sup>Dimensions are measured from the outermost edges of the main frame. Specifications are rounded to the nearest whole number. Conversions are calculated before rounding.

<sup>2</sup>Depth dimensions are measured with extendable handles retracted.

<sup>3</sup>Maximum load capacity is total weight distributed in accordance to basic human anatomy. EMTs must consider the weight of the patient, equipment and accessories when determining the total weight on the product.

<sup>4</sup>Patient restraint set includes two cross-chest restraints and one ankle restraint. Cross-chest restraints can also be used in a chest and lap restraint configuration.

Stryker reserves the right to change specifications without notice.