

# Transfer-Flat

patient transfer

**stryker**<sup>®</sup>



---

## EMS Equipment

3800 E. Centre Ave  
Portage, MI 49002 USA  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 784 4336

[www.ems.stryker.com](http://www.ems.stryker.com)

---

## Standard Features

- Twelve rigid lift handles
- High visibility yellow color
- Compact storage size with carry handle
- 23 oz. vinyl construction
- 2" polyester reinforcement webbing
- Heavy-duty handle reinforcement patches
- Twelve rigid heavy-duty lift handles provide superior versatility and operator comfort
- The Transfer-Flat rolls into compact storage size with convenient carry handle
- Heavy-duty vinyl construction maximizes durability under harsh conditions
- Flexible material facilitates easy transfer for any size patient

---

## Specifications

• <b>Model Number</b>	<b>6083-001-200</b>	
• <b>Flat Size</b>		
Length	108"	(274 cm)
Width	60"	(152 cm)
• <b>Folded Size</b>		
Length	21"	(53 cm)
Width	9"	(23 cm)
Height	9"	(23 cm)
• <b>Weight</b>	10 lbs	(45 kg)
• <b>Maximum Weight Capacity</b>	1600 lbs	(727 kg)

*Stryker reserves the right to change specifications without notice.*