



stryker®

Medical

**Protection and Convenience**  
For You and Your Patient

Secure II®

# Stryker's Secure II® Med/Surg Bed

## The Standard in Patient Care

**Stryker's new siderail design incorporates the necessary changes to ensure compliance with the FDA's dimensional patient entrapment guidance without compromising the exclusive features and benefits that have made us the market leader in patient handling equipment.**

Working in a modern healthcare environment is often a balancing act. Across all units and specialties, a caregiver's ultimate priority is delivering outstanding patient care with momentum and efficiency without sacrificing safety or quality.

That's why Stryker developed the Secure II bed system. Differentiated through innovation, the Secure II is a highly flexible patient care platform designed to provide the ultimate degree of responsive performance.

### Product Highlights:

#### 1. Safe

**Faced with increasing patient weight, age and acuity, clinicians unwittingly take on more risk as they are counted upon to provide superior care despite limited time and resources. Such expectations demand a patient handling solution that not only protects your patients, but also reduces the possibility of compromising your personal safety.**

#### 2. Easy to Use

**Intuitive operation is paramount in today's fast-paced environment where time is critical. A wide range of clinical capabilities coupled with straight-forward, easy to understand controls make the Secure II an exceptionally convenient bed to use.**

#### 3. Proven Design

**Fully customizable through a host of features designed around safety and ease of use, Secure II provides infinitely-adaptable functionality to meet your varying needs. Having a bed system that works for you allows you a greater opportunity to focus on your patient's recovery.**



**Designed for convenience,**  
Engineered for safety.

Open architecture means you have the ability to choose the sleep surface technology most appropriate for your patient, with no loss of bed functionality.





### Fall Prevention

Patient falls are one of the single greatest hazards in today's hospital environment, and preventing them is a daily challenge. Equipment that anticipates caregivers' needs and safety education programs that highlight awareness can help mitigate the risk of patient falls.

The combination of Stryker's Chaperone bed exit system and Fall Prevention Program form the most comprehensive fall prevention solution available.



### Safe

Retractable bed frame eliminates patient's need to reach for belongings, which may increase the risk of falling. Intermediate siderail position provides dual arm support for side ingress and egress. Chaperone® bed exit system senses patient positioning with greater accuracy and reliability.

### Easy to Use

Foot-end nursing controls are easy to access and simple to operate. Integrated scale system provides accurate, repeatable weights with the touch of a button. Centrally-located brake allows uncompromised access with the bed and the siderails in any position. Bed frame is compatible with existing fracture frame equipment and patient helpers.

### Proven Design

Exceptional low height promotes safer patient entry and exit. Standard 85-inch long mattress deck accommodates a 10-inch bed extender with no loss of bed functionality. Low chair position is ideal for patient ambulation. Zoom motorized drive system eliminates pushing.





### Low Chair Position

Secure II's low chair function helps position patients up in bed while maintaining the low bed height necessary to ensure their safety.

### Retractability


Retractable bed frames shorten the bed when the head section is raised. This keeps the patient in a position to access personal belongings stored on bedside tables without the need to reach and risk falling, or call for nursing assistance.

Non-retractable bed frames pull the patient away from bedside tables, which may increase the risk of bed-related falls and calls to nursing staff for help.



### Integrated Scale

Stryker's built-in scale system takes accurate, repeatable weight readings with the patient in any position.



### State-of-the-art Fall Prevention System

Chaperone constantly tracks patient positioning using patented load cell technology. Zone Control® customizes monitoring needs to the patient's level of fall risk.



### Intuitive Siderail Controls

Bed and patient controls located at both head-end siderails and at the foot-end of the bed provide a convenient, three-sided point of care.



### 110-volt Outlet

Conveniently located at the foot-end of the bed, the 110-volt outlet provides an additional source of power to reduce electrical connections to the wall.



### Intermediate Siderail Position

Eliminates the need to reposition the patient up in bed and provides dual arm support for side ingress and egress.



### Bed Extender

Standard 85-inch long mattress deck accommodates a 10-inch bed extender with no loss of bed functionality.



### Centrally-located Brake

Positioned at the middle of the bed, easy-to-engage brakes allow uncompromised access with the bed and the siderails in any position.



### Dual-pedestal Design

Unique bed frame construction eliminates pinch and crush points, and provides supreme cleanliness.



### Zoom® Motorized Transport

Optional Zoom drive system minimizes the effort needed for patient transport while encouraging proper ergonomic posture.

---

**Joint Replacements**

---

**Trauma**

---

**Spine**

---

**Micro Implants**

---

**Orthobiologics**

---

**Instruments**

---

**Interventional Pain**

---

**Navigation**

---

**Endoscopy**

---

**Communications**

---

**Imaging**

---

**Patient Handling Equipment**

---

**EMS Equipment**

6300 South Sprinkle Road  
Kalamazoo, MI 49001 USA  
t: 269 324 2100 f: 269 329 2311  
toll free: 800 787 9537

[www.stryker.com](http://www.stryker.com)

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with <sup>™</sup> designation are trademarks of Stryker.  
Products referenced with <sup>®</sup> designation are registered trademarks of Stryker.

LLB-453-02A  
PC/SCG-15M-0805  
Mkt Lit-114-07282005 Rev A

Copyright © 2005 Stryker  
Printed in USA