

# Symmetry® II

## Bariatric Chair



### **Sleek and Sturdy for Your Bariatric Needs**

Durability and design unite to offer a high-performing Bariatric Chair. Transitional styling and clean lines provide a classic look that is easily integrated in all lounge and patient room environments.

---

**750 lbs dynamic weight capacity**

---

**Removable seat and slip-covered cushions facilitate maintenance, cleaning and enhance durability**

---

**High-density Ultracell foam promotes maximum comfort and resiliency**

---

**Perma-Mesh™ suspension with 2" steel frame provides additional support**

---

**Accentuated lumbar design supports lower back**

---

**Arched, open-arm design with extended arm caps permits easy entry and exit**

---

**Limited lifetime warranty**

# Symmetry® II

## Bariatric Chair

stryker®



---

### Medical

3800 E. Centre Ave.  
Portage, MI 49002  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

[www.stryker.com](http://www.stryker.com)

---

### Standard Features

- Meets flammability standards of CAL 117
- Meets ANSI/BIFMA X5.1-2002 per tests conducted in a certified test facility

---

### Specifications

	Bariatric Chair Part #3500-000-420	
• Overall Width	35"	(89 cm)
• Overall Height	37"	(94 cm)
• Depth		
Overall	25"	(64 cm)
Seat	18"	(46 cm)
• Weight	70 lbs	(32 kg)
• Seat Width	31.5"	(80 cm)
• Back Height From Seat	18.5"	(47 cm)
• Seat Height From Floor	18.5"	(47 cm)
• Arm Height From Floor	26.5"	(67 cm)

*Stryker reserves the right to change specifications without notice.  
Furniture is manufactured exclusively for Stryker Medical.*

---

### Warranty

- Limited lifetime warranty
- Warrantied 750 lbs (340 kg) weight capacity

---

### Optional Features

- Choose from a variety of stains and fabrics
- CAL 133

Copyright © 2006 Stryker  
LSS-457-17D  
PC/SCG-7.5M-0606, 06-SMD-1292  
Mkt Lit-100-02022005 Rev D