

Examination Light



Energy efficient halogen lamp has a rated life of 1,000 hours

Fully-pivoting 360-degree head and spring-loaded arm make positioning easy and accurate

Free-standing portable stand on locking casters offers complete mobility

Specially-designed reflector system ensures a precise light source that renders true colors

Service life backed by two-year warranty

True Light Brings Real Efficiency

Introduce the power of light into your care practices with Stryker's portable examination light. Mobile base and fully-pivoting head allow flexible positioning capabilities to maximize clinician access and visibility. The locking casters and 360-degree spring-loaded articulating arm offer precise positioning without drifting.

Examination Light



Standard Features

- Compact and enclosed articulation system with counterbalance
- Spring-loaded articulating arm
- Sterilizable handle
- 50 Watt halogen bulb

Specifications

• Part #8500-000-001

• Height

Floor to Elbow	69"	(1750 mm)
----------------	-----	-----------

• Length

Elbow joint to center of head	44"	(1117 mm)
-------------------------------	-----	-----------

Castor Base	22"	(565 mm)
-------------	-----	----------

• Width

Castor Base	22"	(558 mm)
-------------	-----	----------

• Weight

	55.8 lbs	(25.3 kg)
--	----------	-----------

• Diameter

Light Head	14"	(355 mm)
------------	-----	----------

• Lamp

Rated Life	1000 hrs
------------	----------

Color Temperature	3978K
-------------------	-------

Color Rendering Index	92
-----------------------	----

• Light Distribution

Luminous Power	33480 fc	(3750 lux)
----------------	----------	------------

Beam Diameter	7" @ 1 m
---------------	----------

• Electrical

Volts	115 V
-------	-------

Frequency	60 HZ
-----------	-------

Power Cord	15'-0" terminating in 3-prong hospital grade plug
------------	---

Switching Power	On/Off
-----------------	--------

Stryker reserves the right to change specifications without notice.

Warranty

- Two year parts and labor

stryker[®]

Medical

3800 E. Centre Ave.
Portage, MI 49002
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com