

LD304

Birthing Bed



Patented Lock-Rite™ foot section glides off for quick and easy breakdown, while a visual indicator gives you confidence it's locked

Length-adjustable footrests provide infinite leg positioning capabilities

Ultra-comfort, 5" thick mattress with integrated, high-memory egg crate topper

Patented Skoocher™ feature moves mom closer to the perineal edge of the bed with the touch of a button

Your Partner in Delivery

Tailor-made for the delivery room, the LD304 birthing bed provides the highest level of comfort and security for mother and newborn. Its enhanced ergonomics, simple operation and advanced mobility also deliver unsurpassed safety and support for caregivers.



Intuitive controls, such as single-button Trendelenburg and Stryker's powered Skoocher feature, afford maximum positioning with minimal effort.



The ergonomic design of Stryker's foot section promotes proper body mechanics to help protect you from back strain and injury.

LD304

Birthing Bed



stryker®

Medical

3800 E. Centre Ave.
 Portage, MI 49002
 t: 269 329 2100 f: 269 329 2311
 toll free: 800 787 9537

www.stryker.com

Standard Features

- Lock-Rite™ slide-off foot section
- Lightweight, self-standing foot section
- Glideaway foot supports with six-way adjustability
- Single-button Trend.
- Three-position siderails
- In-rail motion controls with lockout
- 17.5" (44 cm) low bed height
- 5" (13 cm) Ultra Comfort head and seat mattress
- 3" (8 cm) foot mattress
- 6" (15 cm) Omni Surface™ casters with wheel covers
- Locking-caster steering
- Four-wheel, steel-ring brake system with centrally-located activation pedal
- 33" (84 cm) mattress width with four motor functions (head, bed, foot, Trend.)
- Automatic pelvic tilt
- CPR release
- Manual backup
- Self-storing handgrips
- Plastic fluid basin
- Under-bed night light
- Wood headboard
- Permanent folding three-stage IV pole
- Two year parts and labor warranty

Specifications

• Model Number	4701 (International 4701 model has 230 volt option)	
• Overall Length	93"	(236 cm)
• Overall Width		
Siderails Up	39"	(99 cm)
Siderails Down	37"	(94 cm)
• Lift Capacity		
Bed	500 lbs	(228 kg)
Head Section	250 lbs	(114 kg)
Foot Section	300 lbs	(136 kg)
• Height Range (to litter top)		
High	35.5"	(90 cm)
Low	17.5"	(44 cm)
• Litter Positioning		
Backrest	0° - 70°	
Pelvic Tilt	4°	
Maximum Trend.	8°	
• Patient Surface	33" x 81"	(84 x 206 cm)
• Mattress Thickness		
Head/Seat	5"	(13 cm)
Foot	3"	(8 cm)
• Seat Depth		
Standard	12.5"	(32 cm)
With Skoocher	12.5" - 8.5"	(32 - 22 cm)
• Caster Diameter		
Standard	6"	(15 cm)
Optional	8"	(20 cm)

Electronics

• Agency Approvals	C-US UL (120 VAC), CE (230 VAC)
• Volts	120 VAC (Domestic) / 230 VAC (International)
• Ampere Rating	10A
• Current Leakage	<300 microamperes
• Frequency	60 Hz (Domestic) / 50/60 Hz (International)
• Hospital Grade Plug	Yes
• EMC Compatible	Yes
• Motor Function	Bed, head, foot, Trend., Skoocher (optional)

Stryker reserves the right to change specifications without notice.

Optional Features

- Patented Skoocher™ adjusts seat depth from 12.5" to 8.5" (32 cm to 22 cm)
- Choice of removable or permanent calf supports
- Lumbar pillow and seat pillow
- Mattress overlay
- Basin liners
- Wallsaver™ quick-release connecting kit
- Siderail communications—includes nurse call only or nurse call with speaker, TV, radio, volume, room/read lights
- Smart TV—includes channel up/down, mute
- Pendant control
- Pendant port
- 8" (20 cm) casters
- Labor bar
- Linen for head and foot mattresses
- Extended warranty

Copyright © 2005 Stryker
 LSS-453-03D
 PC/SCG-5M-0206
 Mkt Lit-75-02022005 Rev D