

iBed Awareness™ Option

Medical

3800 E. Centre Ave.
Portage, MI 49002
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com



Monitors brakes, siderail position, bed exit and low bed height, and provides proactive alert if parameters change

Only feature that provides a visual cue when the bed is in a safe state according to the parameters set by the caregiver

Increases nursing efficiency by helping to reduce time associated with routine bed checks and shift changes

Provides head-of-bed lock-out button which prohibits the patient from moving the fowler below 30-degrees

Available with optional audible alarm and ability to connect to nurse call systems

Safe Alert Monitoring

Nurses have more responsibilities today than ever before — Stryker’s iBed Awareness option offers enhanced efficiency that allows them to focus more on care by reducing time spent on routine tasks such as bed checks and shift changes. The iBed Awareness option is an exclusive feature that allows caregivers to select specific parameters that indicate a safe bed for a particular patient. It then locally notifies the caregiver (via a visual and optional audible alarm) if the pre-set parameters change, and the bed is now in an “unsafe” state.



Set the parameters of the bed and press “Bed Status/On” to activate; a green light indicates the bed is “safe”.



When one of the parameters changes, the control panel displays what has been altered and an amber light indicates the bed is no longer “safe”.