

# Stretcher Chairs

Medical



Model 5050

## Increase Patient Throughput

Ideal for surgery centers and small spaces, Stryker's stretcher chairs maximize patient throughput without compromising patient care and outcomes. Choose optimal positioning from horizontal to upright to achieve superior clinical performance and patient comfort.

**Stretcher easily articulates from a flat to a chair position**

**Pneumatic backrest maneuvers 0 to 90 degrees**

**Dual hydraulic jacks and four-wheel, steel-ring brakes optimize stability**

**Foot section can be moved dependently or independently of backrest**

**Eye stretcher chair model features dual-articulation headpiece for eye surgery**



Model 5051

The eye stretcher chair provides all the benefits of the standard stretcher chair plus a dual-articulation, enhanced-clearance headpiece for surgical positioning.

Accent Color Options



# Stretcher Chairs

**stryker**<sup>®</sup>



Model 5050



Model 5051



Model 5051 Supine Position

## Medical

3800 E. Centre Ave.  
Portage, MI 49002  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

www.stryker.com

### Standard Features

- 2" (5 cm) vinyl mattress
- 6" (15 cm) casters
- 24" (61 cm) wide patient surface
- Adjustable height footrest
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Fixed position pushbar (model 5050 only)
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual side activators
- Hands-free Trend.
- Swing-down siderails with comfort armrest
- Thermoformed ABS hood with integrated storage tray and oxygen bottle holder

### Model 5051: All Standard Features Plus

- Dual-articulation, enhanced-clearance headpiece
- Head pads—Ultra Comfort
  - 3" (8 cm), 2" (5 cm) or 1" (2.5 cm) flat
  - 3" (8 cm) or 2" (5 cm) concave

### Specifications

• <b>Model Numbers</b>	<b>5050/5051</b>	
• <b>Overall Length</b>	76"	(193 cm)
• <b>Width</b>		
Overall	30"	(76 cm)
Litter	27.5"	(70 cm)
• <b>Weight Capacity</b>	400 lbs	(182 kg)
• <b>Height Range (to litter top)</b>		
High	33.5"	(85 cm)
Low	22"	(56 cm)
• <b>Litter Positioning</b>		
Backrest	0° - 90°	
Leg	0° - 80°	
Trend./reverse Trend.	±18°	
• <b>Siderails</b>		
Standard	10" x 31"	(25 x 79 cm)
Transfer	9" x 31"	(23 x 79 cm)
• <b>Patient Surface</b>	24" x 74"	(61 x 188 cm)
• <b>Caster Diameter</b>	6"	(15 cm)

*Stryker reserves the right to change specifications without notice.*

### Optional Features

- 2" (5 cm) or 3" (8 cm) Enhanced Comfort Mattress
- 3" (8 cm) Ultra Comfort Mattress
- Concave head pad (model 5051 only)
- Independent foot section
- IV caddy
- IV poles
  - Tethered
  - Standard removable
- Heel stirrups and brackets
- Leg supports and brackets
- Multi-position pushbar (model 5050 only)
- Integrated patient transfer/arm board
- Restraint straps
- Extended warranty available

### Model 5051: All Optional Features Plus

- 3" (8 cm) inflatable head pad
- Adjustable arm board and clamp (surgery accessory rail required)
- Air delivery system/drape support
- Surgery accessory rail (not available with patient transfer siderails)
- Temporal or superior wristrest

Copyright © 2005 Stryker  
LSS-455-04C  
PC/SCG-5M-0206  
Mkt Lit-64-01262005 Rev C

CE