

# Trauma Stretcher

**Full-length X-ray capability**

**Securely-mounted, lateral C-spine  
cassette holder**

**Precisely position slide-out tray for  
supreme imaging results**

**90-degree backrest helps facilitate  
imaging access**



## Stability for the Critically Injured

Stryker's trauma stretcher provides a secure platform for assessing and tending to critically injured and immobile patients. Its innovative design allows for full-length X-ray capability without having to transfer or reposition the patient.



Slide-out cassette tray allows you to capture precise X-ray images without moving the patient.

## Accent Color Options



# Trauma Stretcher

**stryker**<sup>®</sup>



## Medical

3800 E. Centre Ave.  
Portage, MI 49002  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

www.stryker.com

### Standard Features

- 3" (8 cm) mattress
- 8" (20 cm) Omni Surface™ casters
- 26" (66 cm) wide patient surface
- Dual cylinder pneumatic backrest
- Dual mounted collection device hooks
- Dual pedestal hydraulics
- Dual side-mounted foot control hydraulics
- Fibersin pneumatic backrest/stationary foot-end litter
- Fold-down siderails
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual end activators
- Full-length X-ray cassette tray
- Integrated storage tray with oxygen bottle holder
- Roller bumpers
- X-ray centering label
- Constant descent hydraulics
- Retractable fifth-wheel steering
- Wheel covers

### Specifications

|                                       |             |               |
|---------------------------------------|-------------|---------------|
| • <b>Model Number</b>                 | <b>1027</b> |               |
| • <b>Overall Length</b>               | 83"         | (211 cm)      |
| • <b>Overall Width</b>                | 31.5"       | (80 cm)       |
| • <b>Weight Capacity</b>              | 500 lbs     | (228 kg)      |
| • <b>Height Range (to litter top)</b> |             |               |
| High                                  | 38.5"       | (98 cm)       |
| Low                                   | 24.75"      | (63 cm)       |
| • <b>Litter Positioning</b>           |             |               |
| Backrest                              | 0° - 90°    |               |
| Trend./Reverse Trend.                 | ±18°        |               |
| • <b>Patient Surface</b>              | 26" x 76"   | (66 x 193 cm) |
| • <b>Siderails</b>                    | 13" x 55"   | (33 x 140 cm) |
| • <b>Caster Diameter</b>              | 8"          | (20 cm)       |

*Stryker reserves the right to change specifications without notice.*

### Optional Features

- Big Wheel®
- 3" (8 cm) or 4" (10 cm) Enhanced Comfort Mattress
- 4" (10 cm) mattress
- C-spine cassette film holder
- Defibrillator tray/footboard/extender
- Footboard
- Footboard/chart holder
- Backrest X-ray cassette holder
- Heel stirrups and accessory bracket
- IV caddy
- IV poles
  - Permanent folding two- or three-stage
  - Tethered
  - Standard removable
- Patient restraints
- Pop-up push handles
- Serving/instrument tray
- Siderail pads
- Upright oxygen bottle holder
- Integrated patient transfer board
- X-ray tray centering guide
- X-ray grid and holder
- Extended warranty available

Copyright © 2005 Stryker  
LSS-455-08D  
PC/SCG-7.5M-1206  
Mkt Lit-68-01262005 Rev D

CE