

REM


Pressure Redistribution

High-quality mattress replacement system offers varying levels of protection against skin breakdown

Provides maximum patient comfort by reducing pressure and evenly distributing weight

Surface cover is fluid and bacterial resistant


Standard five-year product warranty




REM 3000

Versatile Foam Options


The REM series of sleep surfaces provides patients with maximum body support, comfort and pressure redistribution. Easy to clean and latex-free, REM comes in three models ranging from an economical option for patients not at high risk for skin breakdown to a high-end option that delivers superior pressure redistribution.



REM 1000



REM 2000



1000

2000

3000

Medical

3800 E. Centre Ave.
Portage, MI 49002
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- Supreme wearability
- Improved interface pressures
- Advanced comfort
- Fluid and bacterial resistant cover
- Exceptional cleanability
- Latex free
- Meets flammability standard CAL 117

REM 1000: All Standard Features Plus

- Nylon top with low-slip bottom cover
- A comfortable support surface of high-lofted, convoluted foam
- Head-to-foot reversibility lengthens mattress life

REM 2000: All Standard Features Plus

- Standard Nylon top with low-slip bottom cover (Dartex® top optional)
- Waffle-patterned, three-layer foam provides added pressure reduction
- Edge reinforced by perimeter foam
- A 2.5" (6 cm) urethane base provides proper orthopedic support
- Head-to-foot reversibility lengthens mattress life

REM 3000: All Standard Features Plus

- Standard Dartex top with low-slip bottom cover
- Maximum support and pressure reduction with four highly-resilient foams
- Viscoelastic foam top adjusts to the contours of the body
- Five posture zones in the core foam, including a replaceable center zone

Optional Features

- Custom sizes available upon request (standard lead times do not apply)
- Fire barrier available to meet CAL 129 and Boston IX-11

Specifications

• Model Number	1000	2000/3000
• Overall Length	84" (213 cm)	84" (213 cm)
• Overall Width	35" (89 cm)	35" (89 cm)
• Thickness	6.25" (16 cm)	6" (15 cm)
• Weight Capacity	500 lbs (228 kg)	500 lbs (228 kg)

Stryker reserves the right to change specifications without notice.

Sleep surfaces are manufactured for Stryker Medical.

Dartex® is a trademark of Dartex Coatings, Ltd.

Warranty

- Five-year, non-prorated limited warranty