

Symmetric Aire™

Pressure Redistribution

Promotes wound prevention and early intervention

Internal air reservoirs eliminate external air valves that can cause maintenance issues

Low pre-set valve pressures auto-adjust to patient size and weight

Separate zones for head, torso, leg and heel support individual patient needs



Customized Therapy

Symmetric Aire helps reduce surface pressure that can lead to skin breakdown. This sleep surface features nine dynamic air cylinders in four main anatomical zones that automatically adjust to patient size and weight. Symmetric Aire's reliance on an internal air supply means there are no cords or pumps to get in the way, and no external air valves that can clog and malfunction.

Symmetric Aire™

Pressure Redistribution

stryker®



Medical

3800 E. Centre Ave.
Portage, MI 49002
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

Standard Features

- Internal air reservoir
- Customized therapy
- Sloped heel zone
- Firm foam perimeter
- Sculpted foam topper
- Fluid and bacterial resistant cover
- Latex free

Specifications

• Model Number	2750	
• Overall Length	84"	(213 cm)
• Overall Width	35"	(89 cm)
• Thickness	7"	(18 cm)
• Weight Capacity	500 lbs	(228 kg)

Stryker reserves the right to change specifications without notice.

Sleep surfaces are manufactured for Stryker Medical.

Dartex® is a trademark of Dartex Coatings, Ltd.

Warranty

- Three-year full warranty; five-year prorated warranty

Optional Features

- LVT nylon cover provides extra shear protection; facilitates patient transfers
- Dartex® cover with four-way stretch reduces shear and hammocking; urethane coating increases durability