

# Smart solutions.

## Color Cuff® Reusable Tourniquets

**Designed for durability, strength, and extended use**

Color Cuff reusable tourniquets feature a BarbLock™ retainer which enhances connector security by locking the connector to the hose.

Heavy duty urethane tubing is combined with stainless steel connectors for extended use.

Color Cuff is available in single or dual port configurations with stainless steel Quick Connect or Luer Lock connectors.

**To order the Reusable Color Cuff or additional Stryker tourniquet supplies, contact your Stryker sales representative or call Stryker Surgical Customer Service at 800.253.3210**



### Color Cuff Cylindrical Tourniquets & Accessories

#### Single Bladder

Latex free with 6" tubing and limb protection sleeve

SIZE	COLOR	SINGLE PORT	DUAL PORT
8" x 2.5"	Royal Blue	5922-008-135	5922-008-235
12" x 3"	Gray	5922-012-135	5922-012-235
15" x 3.25"	Green	5922-015-135	5922-015-235
18" x 3"	Red	5922-018-135	5922-018-235
18" x 4"	Red	5922-218-135	5922-218-235
24" X 4"	Yellow	5922-024-135	5922-024-235
30" x 4"	Royal Blue	5922-030-135	5922-030-235
34" x 4"	Purple	5922-034-135	5922-034-235

#### Dual Bladder

Latex free with 6" tubing and limb protection sleeve

18" x 5.5"	Royal Blue	5922-018-145	5922-018-245
24" x 5.5"	Navy Blue	5922-024-145	5922-024-245

\* Luer Lock options available for both Single and Double Bladder.

#### Color Cuff Accessories

	SIZE	CATALOG#
Y Adapter (one each)* <i>Adapts dual port Quick Connector cuff to Quick Connect fill-line</i>	5"	5920-000-304
Y Adapter (one each)* <i>Adapts dual port to Luer Lock to Quick Connect fill-line</i>	5"	5920-000-303
Adapter (one each)* <i>Adapts single port Quick Connect to Luer Lock</i>	5"	5920-000-301
Adapter (set)* <i>Adapts single port Quick Connect to Luer Lock</i>	5"	5920-000-305
Non-Sterile Limb Protector Sleeves <i>(50 sleeves/3 case, 5 individual pouches of 10 each)</i>	6" x 9"	5920-000-104

\*Standard hose color: clear

Available in dual set: one each clear and black hose



---

**Joint Replacements**

---

**Trauma, Extremities, Deformities**

---

**Craniomaxillofacial**

---

**Spine**

---

**Biologics**

---

**Surgical Products**

---

**Neuro & ENT**

---

**Interventional Pain**

---

**Navigation**

---

**Endoscopy**

---

**Communications**

---

**Imaging**

---

**Patient Handling Equipment**

---

**EMS Equipment**

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

**To order the Reusable Color Cuff or additional Stryker tourniquet supplies, contact your Stryker sales representative or call Stryker Surgical Customer Service**  
www.stryker.com | 800.253.3210

4100 East Milham Avenue  
Kalamazoo, MI 49001 USA  
t: 269-323-7700 f: 800-999-3811  
toll free: 800-253-3210

Literature Number:  
9100-001-111 Rev. A  
DDM/PS

Copyright © 2007 Stryker  
Printed in USA