


stryker®

Medical

iBed Awareness™

Smart Bed Monitor



Level 1 Connectivity

Smart Bed Monitoring

Through a smart bed monitoring system, *iBed Awareness*[™] communicates local bed status information to help ensure patient safety and improve overall efficiency.




Set the parameters of the bed and press “iBed” to activate; a green light indicates the bed is monitoring your selected parameters.



A flashing amber light indicates when a parameter has been changed and the control panel displays what specifically has been altered.

***iBed Awareness* is a Stryker exclusive.**
Optional audible alarms and integration into a nurse call system are available.

81% of nurses spend at least 1-2 hours per shift checking beds to proactively ensure patient safety.¹




iBed Awareness

- Allows you to focus more on patient care by reducing time spent on routine tasks such as bed checks and shift changes
- Visually reminds you to reset the Chaperone® Bed Exit System after patient exits bed
- Provides lockout button which prevents the bed from moving below 30 degrees

Three Levels of Connectivity



Level 1: More informed care at the bedside



Level 2: Optimized communication to the nurse station



Level 3: Connecting caregivers throughout the hospital



Joint Replacements

Trauma, Extremities & Deformities

Cranio-maxillofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Pain

Navigation

Endoscopy

Communications

Imaging

Patient Handling Equipment

EMS Equipment

Rehabilitation Services

3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com

REFERENCE

1. Stryker Medical Market Research. March 2007.

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with TM designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

LLB-453-06B
PC/SCG-10M-0507, 07-SMD-1708
Mkt Lit-182-050107 Rev B

Copyright © 2007 Stryker
Printed in USA