


Cordless Driver 2

Stryker, the leader in cordless power


Cordless Driver 2 extends the time-honored tradition of power and versatility in the Sports Medicine arena. With the improved Cordless Driver 2, Stryker continues to offer state-of-the-art battery powered equipment that you have come to expect from Stryker Instruments.

Increased speed and torque are combined with proven cordless reliability to provide unsurpassed power and versatility.

- Powerful
Capable of tunneling, drilling and reaming.



- Versatile
Useful for most orthopaedic procedures. Eliminates the need for multiple handpieces.




- Cordless
Exceptional maneuverability.



Cordless Driver 2 Accessories

Sterilization case

- Modular sterilization containers assist in the handling, storage, and sterilization of your Cordless Driver 2 components



4200-465-000
Sterilization Case

Ordering Information

Cordless Driver 2 Handpiece
4200-000-000 Cordless Driver 2 handpiece
(includes the wire collet 4100-062)

Drilling Attachments

4100-062-000 Wire Collet (.028" - .071" diameter)
4100-110-000 AO Small Drill
4100-125-000 Pin Collet (.078" - .125" diameter)
4100-126-000 Adjustable Pin Collet
4100-131-000 1/4" Keyed Drill Chuck
4100-132-000 5/32" Keyed Drill Chuck
4100-133-000 1/4" Keyless Drill Chuck
4100-134-000 1/8" Keyless Drill Chuck
4100-135-000 Hudson®/Modified Trinkle Drill
4100-160-000 Trinkle Drill
4100-225-000 Pin Collet (.126" - .157" diameter)
4100-600-000 1:1 Bur Attachment
4100-700-000 High Speed Bur Attachment

Reaming Attachments

4100-210-000 AO Large Reamer
4100-231-000 1/4" Keyed Reamer Chuck
4100-235-000 Hudson®/Modified Trinkle Reamer
4100-260-000 Trinkle Reamer
4100-260-000s4 DHS/DCS quick lock

Other Attachments

4100-120-000 T-Latch Attachment
4100-400-000 Sagittal Saw Attachment
4100-310-000 Right Angle Drive - AO Large Reamer
4100-335-000 Right Angle Drive, Modified Trinkle Adapter
4100-355-000 Radiolucent Right Angle Drive
4100-410-000 Radiolucent Adapter



Accessories

4110-120-000 Battery Charger, 120 VOLT
4110-112-000 Battery Adapter for use with 4112 Batteries and System 2000 Charger
4100-420-000 Charger Software Upgrade Kit
(Includes chip and installation tool)
4100-131-132 Universal Drill Chuck Key
4110-412-000 Charger Module for 4112 Battery
4110-422-000 Charger Module for 4212 & 4215 Battery
4112-000-000 Battery Pack
4212-000-000 Small Smart Battery
4215-000-000 Large Smart Battery
4110-426-000 Charger Module for 4126-110, & 4222-110
4222-000-000 Small Aseptic Battery Kit for Cordless Driver 2 (inc. non-sterile battery, transfer shield, battery housing)
4226-000-000 Large Aseptic Battery Kit for Cordless Driver 2 (inc. non-sterile battery, transfer shield, battery housing)
4126-110-000 System 5/Cordless Driver 2 Non-sterile Battery/ Large Cordless Driver 2 Non-sterile Battery/ Small
4222-110-000 Aseptic Sterilization Case (for 4222 battery)
4226-450-000 Aseptic Sterilization Case (for 4226 battery)
4100-451-000 Sterilization Case and Tray
4200-465-000 Plastic Sterilization Case

stryker®

Joint Replacements

Trauma, Extremities & Deformities

Craniomaxillofacial

Spine

Biologics

Surgical Products

Neuro & ENT

Interventional Spine

Navigation

Endoscopy

Communications

Imaging

Patient Care & Handling Equipment

EMS Equipment

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker. Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 1000-621-000 Rev. B
UnDe/PS.

Copyright © 2009 Stryker
Printed in USA

Stryker Instruments
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 800 999 3811
toll free: 800 253 3210

www.stryker.com

stryker®

Surgical

Cordless Driver 2

More Speed
More Torque



Expanding your
Cordless Capabilities.

