

stryker®

Orthopaedics

We Have OR Efficiency  
**Pinned Down**

Introducing  
**PINABALL**  
preloadable pin



# PINABALL preloaded pin

---




The PINABALL preloadable pin can help reduce your OR time during knee replacement surgery. Innovative technology allows you to preload fixation pins into most manufacturers' cutting guides to help lessen the fiddle factor and save time during procedures. The PINABALL device may also help prevent potential OR mishaps,\* such as glove punctures, because pins are loaded prior to hand-to-hand transfer during surgery.

\*The Centers for Disease Control and Prevention has estimated that healthcare workers in hospital settings sustain 384,325 percutaneous injuries involving contaminated sharps annually.

(Federal Registers Occupational Exposure to Bloodborne Pathogens; Needlestick and Other Sharps Injuries; Final Rule. - 66:5317-5325.)

- ▶ Patent pending design uses a spring-loaded ceramic ball at the distal end to hold pin in cutting guide
- ▶ Unique surface-hardening process helps to reduce pin shavings



- ▶ Length: 90 mm (3.54 in)
- ▶ Diameter: 3.2 mm (.125 in)
- ▶ Part# 6003-003-090
- ▶ Pins are packaged sterile (20/box)



