

Maximum Performance

Navigation System II



Maximum Performance

Dependable Results


From Casters to Camera...

The Navigation System II has the performance specifications and accessories you expect from a top of the line surgical navigation system. It features Stryker's own digital camera technology for maximum accuracy and working volume. Dual articulating arms give the surgeon and staff maximum flexibility for positioning in the OR. More features include dual monitors, accessible I/O panel, and powerful computational capabilities. Additionally, the Navigation System II provides the ability to import fluoroscopic and endoscopic images, import scans via a MO drive, and integrate with a microscope.

Stryker Tracking Technology:

- Highly accurate and reliable digital camera.
- Large working volume with virtually no jitter.
- Easy to use wireless control from the sterile field with the use of Stryker's Smart instrumentation.
- Stryker engineered - Stryker manufactured.

Stryker Smart Instruments:




Features:


- Dual Articulating Arms
- Two Large 19" Monitors
- DVD Drive
- Storage Baskets
- Wheel Locks
- Accessible I/O Panel
- MO Drive (optional)
- RFID Capability (optional)

Now with an industrial strength computer:

- Intel® Core™ 2 Duo "Wolfdale" processor
- 4 GB RAM
- 250 GB storage
- 6 times faster than before
- One of the most powerful platforms on the market!




OrthoMap ASM & precision
Knee Navigation Software



OrthoMap Modular
Hip Navigation Software



iNtellect Cranial & ENT
Navigation Software



SpineMap 3D
Navigation Software



Stryker Navigation
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700, f: 800 999 3811
toll free: 800 253 3210

www.stryker.com/navigation

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other affiliated entities own, use or have applied for the following trademarks or service marks: OrthoMap, SpineMap and Stryker. All other trademarks are trademarks of their respective owners or holders.

Literature Number: 9100001154 (Rev. A)
VID/PS 1k 09/09

Copyright © 2009 Stryker
Printed in USA