

Cub® Pediatric Crib

Safety. Access.
Ease of Use.

Responsible pediatric care demands a stretcher specifically designed around safety, access and ease of use. The Cub Pediatric Crib supports safe and efficient care practices while the warm, friendly appearance cultivates a **healing environment for small children and their families.**



Optional protective top equipped with intermediate support position on each side.



Retracting backrest reduces the gap against the siderail while swing-open access doors allow easy patient care.



Lowered siderails provide unobstructed access.

Easily-lowered siderails for **unobstructed access**

Swing-open access doors allow quiet, convenient patient care

Retracting backrest helps eliminate head-end gap

Integrated scale system obtains accurate weights with the child in any position*

Retractable fifth wheel helps improve mobility



Red



Blue



Gray



Green



Purple



Burgundy

Multiple bumper color options allow you to customize the stretcher for identification purposes.

*The scale system allows the reading of a patient's weight in pounds (lb) or kilograms (kg) with a ± 0.5 lb (0.22 kg) precision.

Cub[®] Pediatric Crib

3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com



190 Hydraulic Model



190 Fixed Height Model

Standard Features

- 3" (8 cm) vinyl mattress
- Manually-operated retracting backrest
- Manual foot prop
- Siderail with 9" (23 cm) mandatory stop
- Counterbalance system on siderail for easy raising and lowering
- Oxygen bottle holder
- Fixed foot-end rail

Hydraulic Model:

All Standard Features Plus

- 6" (15 cm) Omni Surface™ casters
- Central-locking brake system
- Dual side-mounted foot control hydraulics
- Uni-lower pedal
- Hands-free Trend./reverse Trend.

Specifications

Model Number	190 Hydraulic CE	190 Fixed Height CE
Overall Length & Width		
No Accessory Supports	65.5" x 37" (166 x 94 cm)	65.5" x 37" (166 x 94 cm)
With Accessory Supports	72.25" x 37.5" (184 x 95 cm)	72.25" x 37.5" (184 x 95 cm)
Weight Capacity	400 lbs (181 kg)	400 lbs (181 kg)
Height Range (to litter top)	39" (86.4 cm)	
High	40" (102 cm)	33" (84 cm)
Low	33" (84 cm)	33" (84 cm)
Litter Positioning		
Manual Activated	0° / 25° / 40° / 55°	0° / 25° / 40° / 55°
Pneumatic Assist	0 - 50°	0 - 50°
Trend./Reverse Trend.	±9.5°	0° / 6° / 12° with manual foot prop
Foot Section Angle	0° / 6° / 12°	0° / 6° / 12°
Patient Surface	29.5" x 57" (75 x 145 cm)	29.5" x 57" (75 x 145 cm)
Caster Diameter	6" (15 cm)	5" (13 cm)

Stryker reserves the right to change specifications without notice.

Warranty

One year parts, labor and travel or two years parts

Optional Features

- 4" (10 cm) Ultra Comfort Mattress
- Integrated scale system (scale system is not CE marked)
- Siderail without 9" (23 cm) mandatory stop
- Access doors at all four corners
- Fully-lowerable foot-end rail - provides unobstructed access around entire unit when all siderails are down
- Accessory supports with roller bumpers (Head, foot or both ends)
- Chart shelf
- Retractable fifth-wheel steering
- IV caddy
- Monitor tray (Not compatible at foot-end when scale option is selected)
- Oxygen bottle retaining ring
- IV poles
 - Permanent folding two- or three-stage
 - Standard removable
- Pneumatic assist backrest
- Retractable protective top
- Retractable protective top with intermediate support position
- Upright oxygen bottle holder
- Siderail pads