


Eye Surgery Stretcher

Maximum Surgical Access.
Simple Functionality.

Stryker's eye surgery stretcher offers the **greatest head-end clearance in the industry**, giving surgeons ample leg room to maneuver while performing delicate eye procedures. The reliable work surface and mobile design combine to increase patient security and surgical throughput.



Flexible tube doubles as air delivery and drape support.



Enhanced clearance and dual-articulating headpiece maximizes surgical access.



Yellow



Gray



Light Blue



Red



Navy



Purple

Multiple bumper color options allow you to customize the stretcher for identification purposes.

Enhanced head-end clearance for **supreme surgeon access**

Ideal for either **superior or temporal** approach

Positioning controls tuck neatly within the headpiece and allow full articulation

Optional **contoured or flat head-piece** comfortably supports patient during procedures

Eye Surgery Stretcher

3800 E. Centre Ave.
Portage, MI 49002 USA
t: 269 329 2100 f: 269 329 2311
toll free: 800 787 9537

www.stryker.com



Standard Features

- 3" (8 cm) Enhanced Comfort Mattress
- 8" (20 cm) Omni Surface™ casters with wheel covers
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Dual side hydraulic controls
- Fold-down siderails
- Four IV receptacles
- Central-locking brakes with foot-end brake & steer pedal
- Hands-free Trendelenburg
- Integrated storage tray with oxygen bottle holder
- Retractable fifth-wheel steering
- Rubber bumpers
- Stationary foot-end litter
- Uni-lower pedal
- Dual-articulation, enhanced clearance headpiece

Specifications

Model Number	1089 CE
Overall Length	90" (229 cm)
Overall Width (siderails up)	31.5" (80 cm)
Weight Capacity	500 lbs. (227 kg)
Height Range (to litter top)	
High	34" (86.4 cm)
Low	22.25" (56.5 cm)
Litter Positioning	
Backrest	0° - 90°
Knee Gatch	0° - 30°
Trend./Reverse Trend.	±18°
Headpiece	
Height Adjustment	±4.5" (11 cm)
Headpiece Tilt	±45°
Patient Surface	26" x 87" (66 x 221 cm)
Siderails	13" x 55" (33 x 140 cm)
Caster Diameter	8" (20 cm)

Stryker reserves the right to change specifications without notice.

Optional Features

- 4" (10 cm) Ultra Comfort Mattress
- Pre/post-op shoulder pads assembly (crank backrest only)
- Air delivery system/drape support
- Crank backrest
- Crank knee gatch
- Defibrillator tray/foot extender
- Temporal or superior wristrest
- Three-sided hydraulic controls
- Footboard/chartholder
- Three-sided brake & steer pedals
- Integrated patient transfer/arm board
- IV caddy
- Restraint straps
- Serving/instrument tray
- Siderail pads
- Surgery accessory rail
- Three-position siderails
- Upright oxygen bottle holder
- Oxygen bottle retaining ring
- Serving/instrument tray holder
- Head pads - Ultra Comfort
 - 3" (8 cm), flat
 - 3" (8 cm), concave
- IV poles
 - Permanent folding two- or three-stage at foot end
 - Removable
- Extended warranty available