

# Stretcher Chairs

Versatile. Simple. Compact.

Ideal for surgery centers and small spaces, Stryker's Stretcher Chairs **maximize patient throughput without compromising patient care and outcomes.** Choose optimal positioning from horizontal to upright to achieve superior clinical performance and patient comfort.



**The Eye Stretcher Chair provides all the benefits of the standard Stretcher Chair plus a dual-articulation, enhanced clearance headpiece for surgical positioning.**

**Multiple bumper color options allow you to customize the stretcher for identification purposes.**



Yellow



Gray



Light Blue



Red



Navy



Purple

Stretcher easily articulates from a **flat to a chair position**

Pneumatic backrest maneuvers **0 to 90 degrees**

**Dual hydraulic jacks** and **four-wheel, steel-ring brakes** optimize stability

Foot section can be moved **dependently or independently** of backrest

Eye Stretcher Chair model features **dual-articulation headpiece** for eye surgery

# Stretcher Chair



3800 E. Centre Ave.  
Portage, MI 49002 USA  
t: 269 329 2100 f: 269 329 2311  
toll free: 800 787 9537

www.stryker.com

## Standard Features

- 2" (5 cm) vinyl mattress
- 6" (15 cm) casters
- 24" (61 cm) wide patient surface
- Adjustable height footrest
- Dual cylinder pneumatic backrest
- Dual pedestal hydraulics
- Four IV receptacles
- Four-wheel, steel-ring brakes with dual-side activators
- Hands-free Trendelenburg
- Swing-down siderails with comfort armrests
- Thermoformed ABS hood with integrated storage tray and oxygen bottle holder
- Multi-position pushbar (model 5050 only)
- Independent foot section

## 5051 Standard Features All Optional Features Plus

- Dual-articulation, enhanced-clearance headpiece
- Head pads - Ultra Comfort
  - 3" (8 cm), 2" (5 cm) or 1" (2.5 cm) flat
  - 3" (8 cm) or 2" (5 cm) concave

## Optional Features

- 2" (5 cm) Enhanced Comfort Mattress
- 3" (8 cm) Enhanced Comfort Mattress
- Concave head pad (model 5051 only)
- IV poles
  - Tethered
  - Standard removable
- Restraint straps
- Extended warranty available

## Specifications

Model Number	5050/5051	CE
Overall Length	76" (193 cm)	
Width		
Overall	30" (76 cm)	
Litter	27.5" (70 cm)	
Weight Capacity	400 lbs (182 kg)	
Height Range (to litter top)		
High	33.5" (85 cm)	
Low	22" (56 cm)	
Litter Positioning		
Backrest	0° - 90°	
Leg	0° - 80°	
Trend./reverse Trend.	±18°	
Siderails	±45°	
Standard	10" x 31" (25 x 79 cm)	
Patient Surface	24" x 74" (61 x 188 cm)	
Caster Diameter	6" (15 cm)	

*Stryker reserves the right to change specifications without notice.*

## 5051 All Optional Features Plus

- Adjustable arm board and clamp (surgery accessory rail required)
- Air delivery system/drape support
- Surgery accessory rail (not available with patient transfer siderails)
- Temporal or superior wristrest