


AlloCraft™ C

Cervical Allograft Spacer




- **Monolithic construct**
- **Chamfered edge**
- **Textured surface**
- **Ported graft**

Contact your local Stryker Spine representative for more information.

AlloCraft™ C

AlloCraft™ C is a cervical allograft spacer designed for Anterior Cervical Discectomy and Fusion (ACDF) procedures. This cervical spacer features a chamfered edge for easy insertion into the disc space. Its monolithic construct has the potential to provide immediate and lasting structural support in the interbody space while the center port allows for the application of a bone graft to help ensure bony fusion.



Part Number	Description
6183-2-006	6mm AlloCraft™ C spacer
6183-2-007	7mm AlloCraft™ C spacer
6183-2-008	8mm AlloCraft™ C spacer
6183-2-009	9mm AlloCraft™ C spacer

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Surgeons must always rely on their own clinical judgment when deciding which treatments and procedures to use with patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The marks bearing the symbol ™ are trademarks of Stryker.
The marks bearing the symbol ® are registered trademarks of Stryker.

Literature Number: LALCSS Rev. 2
GC/GS 200 05/05

Copyright © 2005 Stryker
Printed in USA

325 Corporate Drive
Mahwah, NJ 07430
t: 201 831 5000

www.stryker.com