

AlloCraft™ C-Ring

Cervical Allograft Spacer



- **Monolithic construct**
- **Textured surface**
- **Ported graft**

Contact your local Stryker Spine representative for more information.

AlloCraft™ C-Ring

AlloCraft™ C-Ring is a cervical allograft spacer designed for Anterior Cervical Discectomy and Fusion (ACDF) procedures. This implant features a monolithic construct which has the potential to provide immediate and lasting structural support in the interbody space while the center port allows for application of a bone graft to help initiate bony fusion. The AlloCraft C-Ring is available in a wide range of sizes to accommodate most cervical fusion patients.



Part Number	Description
6183-3-005	5mm AlloCraft™ C-Ring spacer
6183-3-006	6mm AlloCraft™ C-Ring spacer
6183-3-007	7mm AlloCraft™ C-Ring spacer
6183-3-008	8mm AlloCraft™ C-Ring spacer
6183-3-009	9mm AlloCraft™ C-Ring spacer
6183-3-010	10mm AlloCraft™ C-Ring spacer
6183-3-011	11mm AlloCraft™ C-Ring spacer

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Surgeons must always rely on their own clinical judgment when deciding which treatments and procedures to use with patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The marks bearing the symbol ™ are trademarks of Stryker.
The marks bearing the symbol ® are registered trademarks of Stryker.

Literature Number: LALCSS-1 Rev. 3
GC/GS 200 05/05

Copyright © 2005 Stryker
Printed in USA

325 Corporate Drive
Mahwah, NJ 07430
t: 201 831 5000

www.stryker.com