

AlloCraft™ PL

Posterior Lumbar Interbody Fusion Spacer



- **Monolithic construct**
- **Bullet shaped nose**
- **Textured surface**
- **Posterior lordosis**

Contact your local Stryker Spine representative for more information.

AlloCraft™ PL

AlloCraft™ PL is a lumbar spacer designed for use in Posterior Lumbar Interbody Fusion (PLIF) procedures. The AlloCraft™ PL features a bullet shaped nose for easy insertion into the disc space. Its monolithic construct has the potential to provide strength during impaction and following implantation. The implant is fashioned with five degrees of posterior lordosis to allow for a more exact anatomical alignment in the lumbar disc space. Allocraft PL is designed for use with Reliance Total PLIF instrumentation.



Part Number	Description
85008	8mm AlloCraft™ PL spacer
85010	10mm AlloCraft™ PL spacer
85012	12mm AlloCraft™ PL spacer
85014	14mm AlloCraft™ PL spacer
85016	16mm AlloCraft™ PL spacer

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Surgeons must always rely on their own clinical judgment when deciding which treatments and procedures to use with patients. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

The marks bearing the symbol ™ are trademarks of Stryker.
The marks bearing the symbol ® are registered trademarks of Stryker.

Literature Number: LALSS1
GC/GS 200 05/05

Copyright © 2005 Stryker
Printed in USA

325 Corporate Drive
Mahwah, NJ 07430
t: 201 831 5000

www.stryker.com