



## 940 Cast Removal System

### Cast Cutter Handpiece

PART #	DESCRIPTION
940	Cast Cutter

### 940 Cast Cutter Blades

PART #	DESCRIPTION
940-23	<b>Ion Nitrided Blade</b> (replaces standard 2 1/2" blade for use with synthetic or plaster casting material) (non-sterile) 1/pkg.
940-25	<b>Stainless Steel Blade</b> (replaces standard 2 1/2" blade) (non-sterile) 4/pkg.
940-35	<b>Titanium Nitrided Blade</b> (replaces standard 2 1/2" blade for use with synthetic casting material) (non-sterile) 2/pkg.
940-50	<b>Cast Cutter Blade Tool</b> (for blade removal on 940 Cast Cutter) 1/pkg.

### Vacuums

PART #	DESCRIPTION
986	<b>CastVac™ with 8 ft. hose and Mobile Stand</b> (for use with 940, 840, and 848 Cast Cutters)
987	<b>CastVac™ with 8 ft. hose and Wall Mount</b> (for use with 940, 840, and 848 Cast Cutters)

### Vacuum Accessories

PART #	DESCRIPTION
986-500	<b>Replacement CastVac™ Bags</b>
986-848	<b>Adapter Cord and Hose</b> (848 Cast Cutter to 986 CastVac™, 10 ft. cord, 8 ft. hose)
986-840	<b>Adapter Hose</b> (848 Cast Cutter to 986 CastVac™, 8 ft. hose)

### Plaster Dispensers & Cast Room Organizers

PART #	DESCRIPTION
301	<b>Hospital Plaster Dispenser, Stainless Steel</b>  (item 301 also includes the following accessories which may be purchased separately)
300-10-1	<b>Disposable Buckets (non-sterile) 1/pkg.</b>
300-30	<b>Plastic Waste Basket 1/pkg.</b>
300-40	<b>Bucket Liners (non-sterile) 100/pkg.</b>

### General Accessories

PART #	DESCRIPTION
82	<b>Small Cast Spreader</b>
83	<b>Large Cast Spreader</b>
840-27	<b>Cast Cutter Wall Bracket</b>

### Cast Walkers & Walking Heels

PART #	DESCRIPTION
93-100	<b>Standard Adult Walking Heel 24/Pkg.</b>
94-100	<b>Child Size Walking Heel 24/Pkg.</b>
98-100	<b>Medium Size Walking Heel 24/Pkg.</b>
99-100	<b>SofTrac™ Sole Walking Heel 24/Pkg.</b>
96	<b>Toe Caps 12/Pkg.</b>
97	<b>Shower Shields 12/Pkg.</b>

**stryker®**  
**INSTRUMENTS**



## The Stryker 940 Cast Removal System

T H E N E W W O R L D S T A N D A R D

**stryker®**  
**INSTRUMENTS**

For more information, contact your local Stryker Sales Representative.

4100 East Milham Ave  
Kalamazoo, MI 49001  
(800)253-3210  
(616)323-7700

45 Years of Cast Room Experience...

...Combined with Today's Technological Advances

940  
CAST SAW

## 940 Cast Cutter

### Features/Benefits

- **Two Speed Control:**
  - **Standard Speed.** Very quiet, yet powerful enough to remove most casts.
  - **High Speed.** Increased speed for quickly removing even the thickest body jackets and hip spicas.
- **Versatile:** May be used with Vacuum or plugged directly into wall.
- **Quick Connect Vacuum Mount**
- **Impact Resistant Housing**
- **Traditional Ergonomic Design:** Ensures comfort and control while cutting.
- **Integral Vacuum Housing:** Reduces overall size while improving balance and mobility.

- **Efficient DC Motor:** With speed feedback, unparalleled power, less heat build up.

940 Cast Cutter shown actual size.

## Vacuums and Cast Room Accessories

Leading the Industry in Comprehensive Cast Room Experience

## 940 CastVac™

### Features/Benefits

- **Detachable Cannister Lid:** Quick and easy emptying of dust collection filter.
- **Tool Bracket:** conveniently stores small cast spreader, scissors, screwdriver, and blade removal tool on back of stand.
- **Five Wheel Mobile Stand:** Provides stability and easy transportation.
- **Easy to operate Three Way Power Switch**
- **Twelve Foot Cord**



- **Quick Change Blade Mount**



- **Swivel Hose Mount:** Allows for easy movement of the 8ft. Vac Hose.
- **Quiet Vacuum:** Reduces noise level of the entire system.
- **Handpiece Control:** Vacuum automatically turns on when Cast Cutter is activated.
- **986 CastVac™ Specs:** 20" x 20" x 34.5" Tall

## Cutting Accessories

940 Cast blades

