

PainPump[®] 2 BlockAid

Advanced Pain Management



Designed Specifically with the
Anesthesiologist in Mind

Perfect for use in Continuous
Peripheral Nerve Blocks

PainPump[®]2 BlockAid

Advanced Pain Management

**Designed Specifically with the
Anesthesiologist in Mind!**

**Perfect for use in Continuous
Peripheral Nerve Blocks.**

Pump Features:

- Disposable, ambulatory infusion pump with contained fluid reservoir
- Ability to reset flow rates and bolus options
- Reprogrammable options are password protected
- 400mL capacity
- Refillable up to 1000mL
- Separate refill port
- Advanced fluid sensing system
- No audible alarms
- Digital Readout screen



**Provides the sophistication of a resettable, electronic
pump with the convenience of a disposable unit:**

- No hassle with charging/recharging batteries
- No need to track pumps leaving the hospital
- No repair costs for damaged capital items
- No need to clean and sterilize after use

Ordering Information

Re-Programmable Pump
575-100 PainPump[®]2 BlockAid

Indications:

PainPump2 is an electromechanical pump that delivers a controlled amount of medication to the patient for the purpose of managing pain. The pump delivers medication using one or both of the following drug delivery profiles: an hourly flow rate and a bolus PCA (patient controlled anesthetic) dosing option. Medications are infused through intramuscular or subcutaneous routes. PainPump2 is also intended for controlled delivery of local anesthetics in close proximity to nerves for postoperative region anesthesia and pain management. Routes of administration may be intraoperative, perineural, or percutaneous

Contraindications:

- The catheter is not intended for intravenous use
- Infusion of blood and blood products
- Infusion of insulin
- Infusion of critical or life-supporting medication whose stoppage, interruption, overdelivery, or underdelivery would likely cause serious injury or death
- Infusion of any solution that is incompatible with the materials of the PainPump2 system. Consult the pharmaceutical manufacturer's precautions and guidelines and refer to the Technical Specifications in these instructions to ensure the medications to be used will not interact with the materials of the PainPump2 and Infusion Set and possibly cause damage or leakage
- Use in infusion regimens by patients who do not possess the mental, physical, or emotional capability to self-administer their therapy, or who are not under the care of a responsible individual



stryker[®]

Joint Replacements

Trauma

Spine

Micro Implants

Orthobiologics

Instruments

Interventional Pain

Navigation

Endoscopy

Communications

Patient Handling Equipment

EMS Equipment

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with [™] designation are trademarks of Stryker. Products referenced with [®] designation are registered trademarks of Stryker.

Literature Number: **9100-001-003 Rev. None**
UnDe/NAC

Copyright © 2004 Stryker
Printed in USA

4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 800 999 3811
toll free: 800 253 3210

www.stryker.com