

Smart solutions.

## Color Cuff® Sterile Disposable Tourniquets

Surgical

Eliminate the risk of cuff-borne cross contamination with sterile, disposable, Color Cuff cylindrical tourniquets.

Color Cuff combines its low profile, roll resistance and contour compliance to deliver uniform tourniquet pressure. Color Cuff is available in single or dual port configurations with Quick Connect or Luer Lock connectors.



### Color Cuff's patented bladder system offers:

- Improved stability, roll resistance, and contour compliance, delivering uniform tourniquet pressure
- Low profile, stretch resistant bladder
- Velcro® straps provide easy cuff positioning and security
- Cuff ties provide a final level of alignment and safety
- Color coded cuffs that provide easy size identification, increased accuracy, and a simplified inventory

**To order the Disposable Color Cuff or additional Stryker tourniquet supplies, contact your Stryker sales representative or call Stryker Surgical Customer Service at 800.253.3210**

### Color Cuff Cylindrical Tourniquets & Accessories

#### Single Bladder

Latex free with 40" tubing and limb protection sleeve

SIZE	COLOR	SINGLE PORT	DUAL PORT
9" x 1.75"	Light Blue	5921-109-135	5921-109-235
9" x 2.75"	Light Blue	5921-209-135	5921-209-235
12" x 2"	Purple	5921-112-135	5921-112-235
12" x 3"	Purple	5921-212-135	5921-212-235
15" x 2.25"	Pink	5921-115-135	5921-115-235
18" x 3"	Red	5921-018-135	5921-018-235
18" x 4"	Red	5921-218-135	5921-218-235
24" x 4"	Yellow	5921-024-135	5921-024-235
30" x 4"	Royal Blue	5921-030-135	5921-030-235
34" x 4"	Purple	5921-034-135	5921-034-235
44" x 4"	Navy Blue	5921-044-135	5921-044-235

#### Double Bladder

Latex free with 40" tubing and limb protection sleeve

SIZE	COLOR	SINGLE PORT	DUAL PORT
18" x 5.5"	Royal Blue	5921-018-145	5921-018-245
24" x 5.5"	Navy Blue	5921-024-145	5921-024-245

#### Color Cuff Accessories

	SIZE	PART#
Y Adapter (one each)* <i>Adapts dual port Quick Connect cuff to Quick Connect fill-line</i>	5"	5920-000-304
Y Adapter (one each)* <i>Adapts dual port Luer Lock to Quick Connect fill-line</i>	5"	5920-000-303
Adapter (one each)* <i>Adapts single port Quick Connect to Luer Lock</i>	5"	5920-000-301
Adapter (set)* <i>Adapts single port Quick Connect to Luer Lock</i>	5"	5920-000-305

\*Standard hose color: clear  
Available in dual set: one each clear and black hose

\*Luer Lock options available for both Single and Double Bladder.



---

**Joint Replacements**

---

**Trauma, Extremities, Deformities**

---

**Craniomaxillofacial**

---

**Spine**

---

**Biologics**

---

**Surgical Products**

---

**Neuro & ENT**

---

**Interventional Pain**

---

**Navigation**

---

**Endoscopy**

---

**Communications**

---

**Imaging**

---

**Patient Handling Equipment**

---

**EMS Equipment**

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Color Cuff disposable tourniquets are designed for single patient use. Re-use or re-sterilization is not recommended.

**To order the Disposable Color Cuff or additional Stryker tourniquet supplies, contact your Stryker sales representative or call Stryker Surgical Customer Service**

[www.stryker.com](http://www.stryker.com) | 800.253.3210

4100 East Milham Avenue  
Kalamazoo, MI 49001 USA  
t: 269-323-7700 f: 800-999-3811  
toll free: 800-253-3210

Literature Number:  
9100-001-110 Rev None  
DDM/PS

Copyright © 2007 Stryker  
Printed in USA