

stryker®

Surgical

# When Speed and Accuracy Count...

System 6 Heavy Duty  
Sternum Saw

...choose the trusted  
leader in cordless power.



# System 6

## Heavy Duty Sternum Saw

stryker®

Surgical

### Exceptional Cutting Performance



4107-008-000  
Sternum Blade Guard

#### System 6 Sternum Saw

- Increased speed for a more efficient cut
- Cuts 21% faster than the System 5 Sternum Saw
- More power fueled by a new, longer lasting battery
- Smaller, lighter weight handpiece for improved control and comfort



6207-000-000  
Sternum Saw



- See bur & blade guides for a complete list of available cutting accessories

#### Batteries

- Option to use small battery for decreased weight, or larger battery for increased run time

6212-000-000  
System 6 Small Battery

#### Battery Charger

- Quickly charges System 6 batteries
- Modular design allows for 4 different batteries to be charged at the same time
- Interactive screen allows user to assess the condition of the battery and read battery usage
- Auto-reconditions the battery when battery has been damaged by the heat of the autoclave



6110-120-000  
System 6 Battery Charger

6215-000-000  
System 6 Large Battery

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker. Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 9100-001-058 Rev. None UnDe/P.S.

Copyright © 2006 Stryker  
Printed in USA

Stryker Instruments  
4100 East Milham Avenue  
Kalamazoo, MI 49001 USA  
t: 269 323 7700 f: 800 999 3811  
toll free: 800 253 3210

[www.stryker.com](http://www.stryker.com)