

stryker®

Neuro Spine ENT



IMPROVED

Maestro Drill

- Improved Hose
- Improved Air Flow
- Reliable Motor Start-up

access
to better Outcomes

Maestro Drill

Neuro Spine ENT



Improved Hose

- Reduced weight
- Dual elastomer hose provides durability when cleaning
- Will not allow oil to soak through as seen in a silicone hose
- Added compression spring to the hose assembly to prevent exhaust hose from kinking

Improved Air Flow

- Modified muffler to reduce noise

Reliable Motor start-up

- Motor starting feature increases start-up dependability



5400-218-000
Universal Footpedal



5400-206-000
Maestro Irrigation Console

Order Information and Part Numbers

SP/PD Series Attachments

5100-120-952	SD/PD Series Elite MIS Long Curved
5100-120-950	SD/PD Series Elite MIS Long Straight
5100-120-922	SD/PD Series Elite MIS Medium Curved
5100-120-920	SD/PD Series Elite MIS Medium Straight
5100-120-492	SD/PD Series Elite Extra Long Plus Angled
5100-120-490	SD/PD Series Elite Extra Long Plus Straight
5100-120-485	SD/PD Series Elite Super Long Plus Straight
5100-120-484	SD/PD Series Elite Extra Long Straight w/ Bur Guard
5100-120-482	SD/PD Series Elite Extra Long Angled
5100-120-480	SD/PD Series Elite Extra Long Straight
5100-120-472	SD/PD Series Elite Long Angled
5100-120-470	SD/PD Series Elite Long Straight
5100-120-452	SD/PD Series Elite Medium Angled
5100-120-450	SD/PD Series Elite Medium Straight
5100-120-450S1	SD/PD Series XL Elite Medium Straight
5100-120-452S1	SD/PD Series XL Elite Medium Angled
5100-120-470S1	SD/PD Series XL Elite Long Straight
5100-120-480S1	SD/PD Series XL Elite Extra Long Straight

Maestro PD Series Attachments

5400-210-501	Maestro PD Series Adjustable Depth Drill Guide
5400-210-090	Maestro PD Series Extra Long Straight RX
5400-210-080	Maestro PD Series Long Straight R
5400-210-072	Maestro PD Series Angled AM
5400-210-070	Maestro PD Series Straight AM
5400-210-060	Maestro PD Series Perforator Chuck
5400-210-052	Maestro PD Series Angled M
5400-210-050	Maestro PD Series Straight M
5400-210-020	Maestro PD Series Straight O
5400-210-011	Maestro PD Series Right Angle Attachment

Routers

5400-071-059	D-59 Tapered Router
5400-071-058	D-58 Tapered Router
5400-071-057	D-57 Tapered Router

DuraGuards

5400-210-259	Maestro 25mm Steering
5400-210-059	Maestro 25mm Fixed
5400-210-258	Maestro 16mm Steering
5400-210-058	Maestro 16mm Fixed
5400-210-257	Maestro 12mm Steering
5400-210-057	Maestro 12mm Fixed
5400-210-056	Maestro D Attachment

Accessories

5400-279-000	Maestro Sterilization Case
5400-206-001	Maestro Irrigation Console Cable
5400-200-003	Maestro DISS to ARO Adapter
5400-200-002	Maestro DISS to Schrader Adapter
5400-050-001	CORE Irrigation Cassettes (10/pkg)
5400-006-000	Maestro Diffuser (10/pkg)
5400-005-000	Maestro Oil Cartridge (10/pkg)

The information presented in this brochure is intended to demonstrate the breadth of Stryker product offerings. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with TM designation are trademarks of Stryker. Products referenced with [®] designation are registered trademarks of Stryker.

Literature Number: 9100-927-001 Rev. C
UnDe/P.S.
Copyright © 2011 Stryker
Printed in USA

Stryker® Neuro Spine ENT
4100 East Milham Avenue
Kalamazoo, MI 49001 USA
t: 269 323 7700 f: 269 353 3026
toll free: 800 253 3210
www.stryker.com