

Easy Access to a New Segment



Basic Brilliance

Stryker AHTO™ Irrigation System

The AHTO Irrigation Pump is a streamlined irrigation product that is hassle free and extremely reliable for any general surgical case requiring irrigation. This cordless irrigation console is small, lightweight and extremely intuitive. Interactive flow rate and battery level displays establish a user friendly interface with the surgical staff that equates to efficient setup and turnover of the operating room. The simplicity and functionality of the tube set further exemplifies how beneficial this product is to the surgeon, the patient and the surgical staff.



Product Specialties:

AHTO Irrigation Pump

- Maximum Flow Rate of 4 liters/min.
- Up and down toggle buttons for flow control
- Toggle button with battery icon to indicate battery level
- Extended battery life

AHTO Tube Set without Tip

- Delivered Sterile in box of 6
- Cassette that activates the console upon insertion
- Ergonomic hand piece

AHTO Tube Set with Tip

- Delivered Sterile in box of 6
- Disposable tip and tube set
- Cassette that activates the console upon insertion
- Ergonomic hand piece



Features & Benefits:

- Cordless option alleviates clutter in the OR
- Self Priming System
- Displays for Flow Rate & Battery Life
- Instant activation upon insertion of cassette tube set
- Simple installation reduces setup time

System Specifications

AHTO IRRIGATION SYSTEM	SYSTEM SPECIFICATIONS
Dimensions	7.5" h x 6.0" w x 3.75" d
Flow Rates	Low: 2 L/min, Med: 3 L/min, Hi: 4 L/min
Power Requirement	110 V outlet/domestic, 220 V outlet/international

Product Ordering information

ORDERING INFO.	PRODUCTS
250-070-601	AHTO Irrigation Pump
250-070-600	AHTO Tube Set without Tip
250-070-620	AHTO Tube Set with Tip



Stryker Endoscopy
5900 Optical Court
San Jose, CA 95138
t: 1 800 624 4422 f: 1 800 729 2917

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with TM designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

Literature Number: 1000900516 (Rev. B)
VID/FWG 2.5k 12/04