


stryker®

Endoscopy

On The **Inside Tract**

FlexVision™ C-8
Flexible CystoNephroscope



FlexVision™ C-8

On the Inside Tract

Visual Genius

The Stryker enhanced resolution fiberoptic flexible cystonephroscope, the FlexVision™ C-8, provide surgeons an outstanding diagnostic and therapeutic tool to perform their daily urological procedures. The high density fiber bundle provides an enhanced resolution picture not seen in a fiberoptic scope. The 180° deflection both ways allows mapping of the entire bladder and enhanced treatment with flexible instrumentation. The beveled SuperSlim tip and shaft (5.2mm in diameter) ensure easy urethra and prostate navigation leading to greater patient comfort.

Merging the matchless clarity of the FlexVision™ C-8 with the standard-setting Stryker High Definition camera will provide an unrivaled visual image in the office and operating room. Urologists will find the Stryker flexible cystoscope an exceptional and easy-to-use tool by which to perform their most common procedures.



Intuitive Design

The lightweight and ergonomic form of the FlexVision™ C-8 makes ureteroscopies simpler and more comfortable. Combining smooth deflection buttons with a sleek handle offers a firm, dependable grip during procedures.

Superior Deflection

Optimal accessibility is provided in the bladder through primary deflection.



Up



Down



Matchless Clarity

Detailed brightness in the smallest spaces enhances visibility and increases performance during procedures. By joining the forces of the FlexVision™ C-8 with the Stryker Camera System, a surgeon is equipped with the finest tools available to complete complex operations.

FlexVision™ C-8 Features

- » High density fiber bundle for enhanced resolution and picture quality.
- » Beveled SuperSlim tip and shaft (5.2mm diameter) for easy urethra and prostate navigation.
- » Standard flow working channel (2.2mm) for consistent, enhanced flow even with an instrument.
- » 6 o'clock working channel position for optimal instrument positioning and visualization.
- » Thumb control 180° Up/Down responsive deflection for optimal bladder diagnosis and treatment.

Ordering Information

FlexVision™ C-8
503-888-010

The following products are included:

Cleaning Brush
503-888-011

Leak Tester
503-888-005

Port Seal Valve (1 pkg/5)
503-888-001

Shipping Case
1000-000-111

Optional Accessories

Sterilization/Storage Tray
233-032-800

Baskets

Stryker also provides a complete line of Disposable Baskets. Please contact your local Stryker Representative for selection and pricing

Joint Replacements

Trauma

Spine

Micro Implants

Orthobiologics

Instruments

Interventional Pain

Navigation

Endoscopy

Communications

Imaging

Patient Handling Equipment

EMS Equipment

5900 Optical Court
San Jose, CA 95138
t: 800.624.4422 f: 800.729.2917

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label, and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

Literature number: 1000-900-536A
JJ/BL 5k 7/05

Copyright © 2005 Stryker
Printed in USA