

On The Inside Tract



Rigid Cystoscope


On the Inside Tract

The Urology Solution

Stryker is a leader in visualization solutions for minimally invasive procedures. EndoUrologists may now combine the leading 3-chip camera and scope technology with durable, easy to use rigid hardware and instrumentation. Efficient procedures and clear medical images are now packaged into one Urology solution, offered by one company: Stryker.

Superior Performance

Stryker's rigid optical forceps have been engineered for comfort and speed during routine procedures. Standard Cystoscopy sheaths, with a mercier beak and easy-to-read laser marks, are offered in a full range of sizes, each with a durable and ergonomic design. A complete line of instrumentation is available to compliment the sheaths and vary capabilities during all cystoscopic procedures.



Albarran Bridge

Cysto Bridges

Evaluation Bridge
502-880-200

Single Horn Bridge
502-880-201

Double Horn Bridge
502-880-202

Albarran Bridge, Double Horn
502-880-203

Albarran Bridge, Single Horn Bridge
502-880-206



Cysto Bridges




Albarran Bridge with Flexible Grasping Forcep



Clear and Crisp

The goal of Stryker autoclavable endoscopes is simple: provide a clear, crisp image for Urologists. Combined with Stryker's new high definition (HD) camera technology, an excellent picture is now comfortably accessible. First-rate visualization is ensured through Stryker Rigid Optical instrumentation, a line of hardware created for durability and endurance.



Cysto Visual Obturator

Cystoscope Sheath

Flexible Instruments (41cm)

Flexible Grasping Forcep, 5Fr
505-880-105

Flexible Grasping Forcep, 7Fr
505-880-107

Flexible Grasping Forcep, 9Fr
505-880-109

Flexible Biopsy Cup, 5Fr
505-880-115

Flexible Biopsy Cup, 7Fr
505-880-117

Flexible Biopsy Cup, 9Fr
505-880-119

Flexible Scissor, 5Fr
505-880-125

Flexible Scissor, 7Fr
505-880-127

Flexible Scissor, 9Fr
505-880-129



Flexible Biopsy Cup



Flexible Grasping Forcep



Flexible Scissor

Ordering Information

Autoclavable Scopes

4.0mm 0° Autoclavable Cystoscope
502-880-010

4.0mm 12° Autoclavable Cystoscope
502-880-012

4.0mm 30° Autoclavable Cystoscope
502-880-030

4.0mm 70° Autoclavable Cystoscope
502-880-070

Cystoscope Sheaths and Obturators

17Fr Cystoscope Sheath
502-880-217

17Fr Cysto Standard Obturator
502-880-218

19Fr Cysto Visual Obturator
502-881-217

19Fr Cystoscope Sheath
502-880-219

19Fr Cysto Standard Obturator
502-880-220

19Fr Cysto Visual Obturator
502-881-219

21Fr Cystoscope Sheath
502-880-221

21Fr Cysto Standard Obturator
502-880-222

21Fr Cysto Visual Obturator
502-881-221

23Fr Cystoscope Sheath
502-880-223

23Fr Cysto Standard Obturator
502-880-224

23Fr Cysto Visual Obturator
502-881-223

25Fr Cystoscope Sheath
502-880-225

25Fr Cysto Standard Obturator
502-880-226

25Fr Cysto Visual Obturator
502-881-225

Rigid Instruments (Cysto and Resection)

Rigid Optical Biopsy Forcep
502-880-210

Rigid Optical Grasping Forcep
502-880-211

Rigid Optical Cutting Biopsy Forcep
502-880-212

Rigid Optical Crushing Forcep
502-880-213

Rigid Optical Biopsy Forcep 70°
502-880-214

Rigid Optical Biopsy Forcep 30°
502-880-215

Rigid Optical Scissor Forcep
502-880-216

Adapters

Bladder Syringe Adapter, Cysto Sheath
502-880-204

Ellik Adapter, Cysto Sheath
502-880-601

Port Seal (10 per box)
502-880-205

stryker[®]

Joint Replacements

Trauma

Spine

Micro Implants

Orthobiologics

Instruments

Interventional Pain

Navigation

Endoscopy

Communications

Imaging

Patient Handling Equipment

EMS Equipment

5900 Optical Court
San Jose, CA 95138
t: 800.624.4422 f: 800.729.2917

www.stryker.com

The information presented in this brochure is intended to demonstrate a Stryker product. Always refer to the package insert, product label, and/or user instructions before using any Stryker product. Products may not be available in all markets. Product availability is subject to the regulatory or medical practices that govern individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Products referenced with ™ designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

Literature number: 1000-900-512C
JJ/BL 5k 5/05

Copyright © 2005 Stryker
Printed in USA