

stryker®

Endoscopy

DVD Quality

For Today's O.R.

Premier DVD Recorder



DVD Quality For Today's O.R.

Stryker Premier DVD Recorder

Boost the quality and resolution of your operating-room video with the Stryker Premier DVD Recorder. An affordable alternative to operating-room VCRs, the Premier DVD Recorder burns real-time video for playback on commercial DVD players or personal computers equipped with DVD-ROM. The Premier DVD Recorder features simple, straightforward controls that cut the learning curve of upgrading to new technology, plus it is compatible with any surgical specialty. So whether you're recording video for patient viewing or for teaching purposes, the Stryker Premier DVD Recorder delivers the quality and convenience of DVD at a price you can afford.

Product Features

- Burns moving video in real time
- Records in DVD-video mode for playback on commercial DVD players and home computers*
- Records directly from the camera's DVI input and routes to any DVI monitor
- Provides time-lapse and daily recording options
- Records to any DVD format (DVD±R, DVD±RW)



* DVD-ROM required

Ordering Information

240-020-505 Stryker Premier DVD Recorder

Available Accessories

240-050-811 (10 Pack) Sony DVDs, +RW

240-050-895 (10 Pack) Sony DVDs, +R

240-050-894 (10 Pack) Sony DVDs, -R

System Specifications

Dimensions	8.35" w × 13.87" d × 3.47" h
Weight	10.60 lbs
Power-supply input	100 – 240V AC
Power consumption	1.00 – 0.45 A
Video signals	DVI, line and Y/C
Frequency	50 – 60 Hz
Operating temperature	41 – 95° F
Operating humidity	20 – 80% rh

Products referenced with ™ designation are trademarks of Stryker.
Products referenced with ® designation are registered trademarks of Stryker.

Literature number: 1000-900-550A
JJ/BL 5k 02/06

Copyright © 2006 Stryker
Printed in USA

stryker®

5900 Optical Court
San Jose, CA 95138
t: 800.624.4422 f: 800.729.2917

www.stryker.com