

# Prime Connect™



## Smart stretcher. Confident care.

Empower your team with Prime Connect—Stryker's first smart, connected stretcher. With both visual and audible bed exit alarms, Prime Connect can communicate with the nurse call dome light, a Vocera Badge, or the nurse's station through Vision to help enable fall protocol compliance from anywhere in your facility.

### Enhanced local monitoring

Prime Connect makes it easy to acknowledge bed exit alarms in loud, fast-paced environments with the ability to illuminate both the corridor light and footboard display when activated.

### Central monitoring

Our Vision dashboard displays what you need to know in near real-time, including an at-a-glance view of the stretcher's location, fall-risk status, bed exit status and more.

### Mobile monitoring

For on-the-go monitoring and rapid response times, Prime Connect can send notifications to the right care team members through Vocera Engage.\*

### Advanced mobility options

Big Wheel reduces start-up push force and steering efforts, while Zoom motorized drive nearly eliminates manual pushing.

### Recovery chair position

Helps eliminate the need for unnecessary transfers by providing a comfortable seated position.

### Lift Assist® backrest

Uses the patient's weight to raise the head of the stretcher while the drop seat design helps reduce boosting. Also available as a fully electric option.

### Patient controls

Allows for fewer interruptions for you. The patient has control of their own comfort with the staggered control panels on the Glideaway siderails.

### Electric-powered lift pedal\*\*

Eliminates manual pumping, allowing for quick patient access.

### IsoFlex™ SE gel support surface<sup>1</sup>

Offers pressure redistribution and comfort as a component to aid in the treatment and prevention of pressure injuries.

### Power-washable design\*\*

With 700 lb safe working load.



# Prime Connect

## Standard features

- Nurse call integration
- Wireless capabilities
- Enhanced footboard
- Bed exit
- Integrated scale
- Less than 23-inch low height (Models 1105 and 1115)
- Litter-mounted oxygen bottle holder
- Center-locking brake system with dual end activators
- Dual, side-mounted foot control hydraulics with uni-lower pedal
- Dual end brake/steer controls
- Glideaway siderails
- Four integrated IV receptacles
- Trendelenburg/reverse Trendelenburg positioning
- 8-inch Omni surface Casters with integrated wheel covers
- Pneumatic backrest/stationary foot
- Thermoformed ABS base
- Integrated utility tray
- Retractable cord
- Pop-up push handles
- Electric powered knee gatch
- Lift Assist backrest design with drop seat
- Three-sided hydraulic controls
- Dual-ended siderail release

## Specifications

	30-inch litter
<b>Model number</b>	
5th Wheel Mobility	1105
Big Wheel Advanced Mobility	1115
Zoom Motorized Drive System	1125
<b>Overall length</b>	
1105 and 1115	85" (215.9 cm)
1125	86" (218.4 cm)
<b>Overall width</b>	
Siderails up	38" (96.5 cm)
Siderails down	30.5" (77.5 cm)
<b>Patient surface</b>	
30" x 75.25" (76 x 191 cm)	
<b>Weight capacity</b>	
700 lb (318 kg)	
<b>Height range (1105, 1115)</b>	
High	34" (86.4 cm)
Low	20.75" (57.7 cm)
<b>Height range (1125)</b>	
High	34" (86.4 cm)
Low	23" (58.4 cm)
<b>Litter positioning</b>	
Backrest	0° – 90°
Knee gatch	0 – 5.5" (13.9 cm)
Trend./Reverse Trend. (1105, 1115)	+17/-17 (+/- 3)
Trend./Reverse Trend. (1125)	+17/-16 (+/- 3)
<b>Siderail length</b>	
14" x 57.75" (36 x 147 cm)	
<b>Caster diameter</b>	
8" (20 cm)	
<b>EMC</b>	
IEC 60601-1-2	

## Optional features

- IsoFlex SE support surface
- ComfortGel SE support surface
- UltraComfort SE support surface
- Permanent two-stage IV pole
- Serving tray
- Defibrillator tray
- Defibrillator tray/foot extender/charting surface
- Chartholder
- Serving tray holder/footboard
- Siderail pads
- Upright oxygen bottle holder
- Restraint straps

## Warranty

2 year parts, labor, travel

## Electric functionality

- Electric patient controls
- Foot end operator controls
- Electric lift
- Bed exit
- 120V or 240V available

### Reference:

1. Recommended to be implemented in combination with clinical evaluation of risk factors and skin assessments made by health care professional.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: Big Wheel, ComfortGel, Glideaway, IsoFlex, Lift Assist, Prime Connect, Prime Series, Stryker, Vision, Vocera, Vocera Engage, Zoom. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

\*Vocera Engage contains Engage Medical Device Alarm Notification (EMDAN), FDA 510(k)-cleared middleware, to deliver secondary alarm notifications.

\*\*Models 1105 and 1115 only