

Unity[®] Recliner Series

Refresh your room with our furniture line.



Modern, durable and comfortable recliners are designed to be customizable to meet the needs of each department. Clinical and Treatment models were developed to meet the needs of both caregivers and patients to help aid in the recovery process and allow patients to get into a recliner sooner.

Key features

- Wall saver recline designed to help maximize space in compact patient rooms.
- Clinical and Treatment models feature transfer arm options to allow for caregiver access to help better facilitate safe patient handling, transfers, and early mobility. Armrest rotates up-and-back to provide full clearance when the patient is in the full recline position.
- Performance cushion material creates a supportive and comfortable seat, allowing for extended time periods to be spent in the recliner.
- Extended wave-shaped arms promote independence when transitioning between sitting and standing positions.
- Large dual-tread casters with steer lock helps ensure safety and aids in maneuverability in the room and throughout the facility.
- Full 360° clean out with gap between seat and back.



Specifications

Specifications	Value	Clinical narrow Treatment narrow	Clinical wide Treatment wide
Model number	3800-000-101	3800-000-201 3800-000-301	3800-000-202 3800-000-302
Overall width	31" (78.7 cm)	31" (78.7 cm)	35" (88.9 cm)
Overall height	46" (116.8 cm)	46" (116.8 cm)	47" (119.4 cm)
Depth			
Overall	43" (109.2 cm)	43" (109.2 cm)	43" (109.2 cm)
Seat	19.5" (49.53 cm)	19.5" (49.53 cm)	19.5" (49.53 cm)
Trendelenburg depth		74" (188.0 cm)	76.5" (194.3 cm)
Weight (without accessories)	165 lb (74.8 kg)	165 lb (74.8 kg)	175 lb (79.4 kg)
Seat width	22.5" (57.2 cm)	22.5" (57.2 cm)	27.5" (69.9 cm)
Back height from seat	26.75" (67.9 cm)	26.75" (67.9 cm)	27.25" (69.2 cm)
Seat height from floor	21.5" (54.6 cm)	21.5" (54.6 cm)	21.5" (54.6 cm)
Arm height from floor (highest point)	30" (76.2 cm)	30" (76.2 cm)	30" (76.2 cm)
Safe working load	500 lb (226 kg)	500 lb (226 kg)	500 lb (226 kg)

Model descriptions

Value

Narrow seat width with fixed armrests.

Clinical

Foley bag hooks on both sides and 1 transfer arm on the left or right side.

Option to add an additional transfer arm so that both sides open for transfers.

Offers both a narrow and wide seat width.

Treatment

Foley bag hooks on both sides, 2 hook IV pole mounted on a fixed arm, with a transfer arm on the opposite side. Offers both a narrow and wide seat width.

Warranty

Limited lifetime warranty

Warrantied 500 lb (226 kg) weight capacity

Ratings/Awards/Certifications

Fire rating

- All fabric and foam are tested and certified, meeting California Technical Bulletin 117 (CAL 117) standards
- Meets or exceeds all applicable BIFMA standards for safety and durability

Standard features

- Powder coated steel base
- Wall saver recline allows chair to be fully functional only 8" from wall
- Dual-sided controls
- Transfer arm on both Clinical and Treatment recliner models
- Steer lock caster
- Push handle
- Trendelenburg position on both Clinical and Treatment recliner models
- 4" casters
- 2 Foley bag hooks on both Clinical and Treatment recliner models
- Clearance for Stryker's Tru-Fit® OBT
- 360° cleanout access around seat cushion for easy housekeeping and infection control practices
- Clean and disinfect with chlorinated bleach solutions, isopropyl alcohol, quaternary solutions and Stryker's SideKick® solutions

Optional features and accessories

- Central locking brake
- Additional transfer arm available on Clinical recliner models
- Fabric grades A-C
- Foot tray (500 lb SWL)
- Rear-mounted IV pole on Treatment recliner models
- Head flap
- Head pillow
- Bolster cushion for seat width adjustability