

stryker

LIFEPAK® 15

monitor/defibrillator



Accessories guide

LIFEPAK 15 monitor/defibrillator

Power options



Mobile battery charger
Includes AC and DC power cords, mounting bracket and operating instructions.
11577-000011



Station battery charger
Includes AC and DC power cords, and operating instructions.
11577-000004



Lithium-ion battery 5.7 Ah
5.7 amp hour, 11.1 volt, rechargeable, with fuel gauge.
21330-001176



AC power adapter
Includes Right Angle Cable (AC power cord not included).
11140-000098

DC power adapter
Includes DC cable and Right Angle Cable.
11140-000074

AC power cord
11140-000015
Power attachment kit
11577-000019



Extension cable
For AC/DC Power Adapter.
11140-000080




Replacement right angle power cable
For AC/DC Power Adapter.
11140-000081



Replacement DC input cable
For DC Power Adapter.
11140-000084

ECG monitoring accessories




12-lead ECG cable trunk cable with 4-wire limb leads

11111-000018 (5ft)
11111-000020 (8ft)



12-lead ECG cable 6-wire precordial attachment

11111-000022




4-wire cable comb
21300-008054 (10/pack)


6-wire cable comb
21330-008055 (10/pack)



3-wire ECG cable
11110-000029



5-wire ECG cable
Right-angle connector, 4-wire limb plus one chest lead, labeled "V1" on the LIFEPAK 15 monitor reports.
11110-000066



Box of strip chart recorder paper
100mm x 22m
11240-00003 2 (2 rolls per box)



ECG electrodes, high adhesion

31115796 (200 3/packs, 600 total)
31115788 (120 5/packs, 600 total)

LIFEPAK 15 monitor/defibrillator

Therapy delivery accessories

Hard paddles and electrode gel



Standard hard paddles
11130-000061 (1 pair)



Pediatric paddle, external
(Two required); slips onto standard adult hard paddle.
11133-000007 (1 adapter)



SIGNAGEL® electrode gel
For use with hard paddles. Highly conductive, multipurpose electrolyte meets all the standards of the ideal saline electrode gel. Recommended for ECG, defibrillation, biofeedback and EMG.
21300-005847 (8.5 oz)



Sterilizable internal defibrillation paddles (U.S. only)
For use with LIFEPAK 15 and 20e defibrillator/monitors.

1 inch size
11131-000044 (1 pair)
1.6 inch size
11131-000047 (1 pair)
2.3 inch size
11131-000045 (1 pair)
3 inch size
11131-000046 (1 pair)



Internal paddle handles with discharge control (OUS only)
(For use with the internal paddles adapter cable)
11131-000001 (1 pair)



Internal paddles adapter cable (OUS only)
(For use with internal paddle handles)
11998-000326



Internal paddles (requires internal paddle handles and internal paddles adapter cable) (OUS only)

1 inch size
11131-000010 (1 pair, 6.25 inch shaft)
1.5 inch size
11131-000011 (1 pair, 6 inch shaft)
11131-000021 (1 pair, 9 inch shaft)
11131-000024 (1 pair, 14 inch shaft)
2 inch size
11131-000012 (1 pair, 5.75 inch shaft)
11131-000022 (1 pair, 8.75 inch shaft)
2.5 inch size
11131-000013 (1 pair, 5.5 inch shaft)
11131-000019 (1 pair, 8.5 inch shaft)
3.5 inch size
11131-000014 (1 pair, 5 inch shaft)
11131-000023 (1 pair, 8 inch shaft)

Therapy delivery accessories

EDGE System™ electrodes for pacing/defibrillation/ECG with QUIK-COMBO® connector

18-month minimum shelf life remaining at time of shipment from Stryker except where noted.



EDGE System electrodes with QUIK-COMBO connector
24-inch leadwire length
11996-000091



EDGE System RTS (radiotransparent) electrodes with QUIK-COMBO connector
24-inch leadwire length
11996-000090



EDGE System electrodes with QUIK-COMBO connector and REDI-PAK® preconnect system
24-inch leadwire length
11996-000017



Pediatric EDGE System RTS electrodes with QUIK-COMBO connector

For use only with manual monitor/defibrillators; 12-month minimum shelf life at time of shipment 24-inch leadwire length.
11996-000093



QUIK-COMBO therapy cable

With convenient TRUE-LOCK® cable connector. Length is approximately 8 ft. For use with LIFEPAK 15 monitor/defibrillator.
11113-000004 (8 ft)

LIFEPAK 15 monitor/defibrillator

NIBP monitoring accessories

NIBP hoses



NIBP tubing, straight

21300-008159 (6 ft)
21300-008147 (9 ft)
21300-008146 (12 ft)



NIBP tubing, coiled

21300-008148 (2-9 ft)

NIBP cuffs



Reusable cuff

X-large adult 35 - 44 cm 11160-000019	Pediatric 13 - 20 cm 11160-000013
Large adult 32 - 42 cm 11160-000017	Infant 8 - 14 cm 11160-000011
Adult 2 26 - 35 cm 11160-000015	



Single patient use cuff

X-large adult 35 - 44 cm 11160-000020	Pediatric 13 - 20 cm 11160-000014
Large adult 32 - 42 cm 11160-000018	Infant 8 - 14 cm 11160-000012
Adult 26 - 35 cm 11160-000016	

Pulse oximetry monitoring accessories

Masimo SET® RC patient cables



RC patient cable
For use with M-LNCS and Rainbow patient sensors.
11171-000037 (4ft)
11171-000038 (12ft)




RC EMS patient cable
11171-000082 (4 ft)


Masimo SET RC patient cable compatible SpO₂ sensors



M-LNCS reusable sensor
11171-000046 (Ad)
11171-000047 (Ped)




M-LNCS adhesive sensors (20/box)
11171-000039 (Ad)
11171-000040 (Ped)



M-LNCS adhesive sensors (20/box)
11171-000042 (Neo/Ad)
11171-000041 (Inf)

Masimo SET RC patient cable compatible Rainbow® SpO₂, SpCO, SpMet sensors



Rainbow reusable sensor
11171-000049 (Ad)
11171-000050 (Ped)



Rainbow adhesive sensors (10/box)
11996-000342 (Inf)
11996-000341 (Neo/Ad)



Rainbow adhesive sensors (10/box)
11996-000339 (Ad)
11996-000340 (Ped)

LIFEPAK 15 monitor/defibrillator

Pulse oximetry monitoring accessories

Masimo SET LNC patient cables**Red LNC patient cable**

For use with LNCS patient sensors.

11996-000323 (4 ft)

11996-000324 (10 ft)


Masimo SET LNC patient cable compatible SpO₂ sensors**LNCS[®] reusable sensor**

11171-000017 (Ad)

11171-000018 (Ped)

**LNCS reusable Soft sensor**

11171-000052 (Ad)

**LNCS adhesive sensors (20/box)**

11171-000019 (Ad)

11171-000020 (Ped)

**LNCS adhesive sensors (20/box)**

11171-000029 (Neo/Pt)

11171-000028 (Neo/Ad)

11171-000031 (Inf)

Pulse oximetry monitoring accessories

Direct connect SpO₂ only patient sensors



Adult reusable direct connect sensor
11996-000331 (3 ft)



Pediatric reusable direct connect sensor
11996-000333 (3 ft)




Adult reusable soft direct connect sensor
11171-000053 (8 ft)

Direct connect Rainbow SpO₂, SpCO, SpMet patient sensors




Adult Rainbow direct connect reusable sensor
11996-000335 (3 ft)
11171-000032 (8 ft)



Pediatric Rainbow direct connect reusable sensor
11171-000033 (8 ft)

Additional Masimo accessories



Reusable ambient light shield
11171-000054 (5/bag)



Disposable ambient light shield
11171-000055 (10/bag)

LIFEPAK 15 monitor/defibrillator

Pulse oximetry monitoring accessories

Masimo to Nellcor™ adapter



Red MNC cable
 Connects LIFEPAK 15 to
 Nellcor patient sensor
11996-000365 (4 ft)
11996-000366 (10 ft)

Masimo to Nellcor adapter compatible SpO₂ sensors



Nellcor reusable SpO₂ sensor, adult, DS100A
11996-000060 (Adult - Ref# DS100A)



Nellcor reusable multisite SpO₂ sensor, D-YS
11996-000106 (≥ 1 kg - Ref# D-YS)



Nellcor reusable two-piece SpO₂ sensors, OXI-A/N, OXI-P/I
 Includes 50 disposable adhesive sensors.
11996-000061 (Ad/Neo)
11996-000062 (Ped/Inf)

Nellcor single-patient use adhesive sensor wraps (100/pk)

Not available in Canada.
11996-000048 (Ad/Neo - Ref# ADH A/N)
11996-000049 (Ped/Inf - Ref# ADH P/I)

Nellcor single-patient use foam wrap

Not available in Canada.
11996-000110 (Ad/Neo - Ref# FOAM A/N)
11996-000108 (Ped/Inf - Ref# FOAM P/I)



Nellcor DEC-4 compatible SpO₂ adapter cable
11110-000042 (4 ft)

Nellcor DEC-8 compatible SpO₂ adapter cable
11110-000176 (8 ft)

Temperature monitoring

Temperature cables



Temperature adapter cable

11140-000079 (10 ft)
11140-000078 (5 ft)

LIFEPAK 15 monitor/defibrillator

End-Tidal CO₂ (EtCO₂) monitoring accessories

Microstream™ Advance short-term filter lines

Single patient use



Orange connectors designate short-term use. For use with Microstream capnography monitors and multiparameter monitors equipped with Microstream technology.



Non-intubated filter lines:

- Adult with O₂
- MVAO** (25/pack, 200 cm)
- MVAO100U** (100/pack, 200 cm)
- MVAOL** (25/pack, 400 cm)
- Adult without O₂
- MVA** (25/pk, 200 cm)
- MVA100U** (100/pack, 200 cm)
- Pediatric with O₂
- MVPO** (25/pack, 200 cm)
- Pediatric without O₂
- MVP** (25/pack, 200 cm)



Intubated filter lines:


- Adult/Pediatric
- MVAI** (25/pack, 200 cm)
- MVAIL** (25/pack, 400 cm)
- MVAI100U** (100/pack, 200 cm)

Microstream™ Advance long-term filter lines

Single patient use



Yellow connectors designate long-term use. For use with Microstream capnography monitors and multiparameter monitors equipped with Microstream technology.



Intubated filter lines:

- Adult/Pediatric
- MVAIH** (25/pack, 200 cm)
- Infant/Neonatal
- MVIIH** (25/pack, 200cm)

Cases and mounting options



Standard carrying case
Includes right pouch and left pouch, for use with LIFEPAK 15 monitor/defibrillator.
11577-000002

Top pouch
Storage for sensors and electrodes; insert in place of standard paddles.
11220-000028

Shoulder strap
For use with LIFEPAK 15 monitor/defibrillator.
11577-000001

Back pouch
Ideal for additional accessory storage. For use with LIFEPAK 15 monitor/defibrillator.
11260-000039



Bed connector
For use in hospital only.
11996-000374

Communication accessories



LIFEPAK monitor to PC cable
For connecting LIFEPAK 12 or LIFEPAK 15 monitor/defibrillator to PC.
11996-000369 (USB)



Titan III wireless gateway
For transmitting data from LIFEPAK 15 to the LIFENET[®] system or CODE-STAT[™] data review software. Requires existing wireless network.
21996-000109 (wireless gateway)
11996-000479 (wireless/cellular/audio gateway) - For use with Stryker Verizon data plan
11996-000484 (wireless/cellular/audio gateway) - For use with Stryker AT&T data plan
11996-000480 (wireless/cellular/audio gateway) - For use with customer data plan



4G gateway
For transmitting data from LIFEPAK 15 to the LIFENET system or CODE-STAT data review software. Requires data plan.
11996-000476 (AT&T gateway)
11996-000474 (Verizon gateway)
11996-000471 (Verizon gateway - for use with Stryker data plan)
11996-000475 (AT&T 4G Kore modem/MultiTech cellular only)

Stryker Data Plans
99428-000305 Verizon data plan 1-year
99428-000304 AT&T data plan 1-year

LIFEPAK 15 monitor/defibrillator

Testers and training materials




12-lead patient simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillation, tachycardias, and bradycardias, as well as ST segment and T wave abnormalities. For use with LIFEPAK devices with a 12-lead ECG feature.
11996-000311 (12-Lead)

3-lead patient simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe, interactive training. Select from 17 ECG rhythms including: fibrillation, tachycardias and bradycardias.
11996-000310



Defibrillator checker


Tests integrity of energy delivery through standard hard paddles. Neon light indicates energy has been delivered.
11998-000060



Test load

Use to perform therapy cable performance checks. Connects to the QUIK-COMBO therapy cable on the defibrillator.
21330-001365 (English only)
21330-001367 (English, French for Canada)

Training tools



Operating instructions: LIFEPAK 15 monitor/defibrillator
 Also available as free online download at <https://ifu.stryker.com>
26500-002408

For further information, please contact Stryker at 800 442 1142 (U.S.), 800 668 8323 (Canada) or visit our website at stryker.com

Emergency Care

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CODE-STAT, EDGE System, LIFEPAK, REDI-PAK, QUIK-COMBO, Stryker, TRUE-LOCK, Masimo, the Radical logo, Rainbow and SET are registered trademarks of Masimo Corporation. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

US-CA 04/2023
 Copyright © 2023 Stryker

Distributed by:

Physio-Control, Inc.
 11811 Willows Road NE
 Redmond, WA, 98052 U.S.A.
 Toll free 800 442 1142
stryker.com

Distributed by:

Stryker Canada
 2 Medicorum Place
 Waterdown, Ontario
 L8B 1W2
 Canada
 Toll free 800 668 8323