

LIFEPAK[®] 20e

defibrillator/monitor with
CodeManagement Module[®]

Genuine accessories

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module





Genuine accessories from **Stryker**

Only Stryker accessories and disposables have been thoroughly tested with LIFEPAK products. Stick with the name you trust and the company that stands behind it. To order, contact your Stryker representative.

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

Power options



AC Power Cord

Powers LIFEPAK 20e defibrillator/
monitor and recharges the
internal battery.

11140-000015



Rechargeable Internal Battery

For use as AC power backup.
11141-000112 (LIFEPAK 20e Lithium-ion)



Lithium-ion Battery for CodeManagement Module

Rechargeable battery for the CodeManagement Module;
for use as AC power backup.
11141-000162 (Lithium-ion)

ECG monitoring accessories



5-Wire ECG Cable
11110-000066



3-Wire ECG Cable
11110-000029



**530-series Medi-Trace®
ECG electrode**

Comes in minimum quantity 600 ECG electrodes.

31115796 (200 packs of 3 ECG electrodes)
31115788 (120 packs of 5 ECG electrodes)



50mm x 22mm ECG Printer Paper
11240-000031 (3 rolls/box)



**CodeManagement Module for use
with LIFEPAK 20e Defibrillator/
Monitor - Wireless**

Connects to the LIFEPAK 20e defibrillator to add wireless connectivity.

11150-000018

**CodeManagement Module for use
with LIFEPAK 20e Defibrillator/
Monitor - Wireless & Capnography**

Connects to the LIFEPAK 20e defibrillator to add wireless connectivity and capnography monitoring.

11150-000019

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

Therapy delivery accessories

Hard paddles and electrode gel



SIGNAGEL® Electrode Gel

For use with hard paddles. Highly conductive, multi-purpose electrolyte meets all the standards of the ideal saline electrode gel.

21300-005847 (8.5 oz)

Standard Adult Detachable Hard Paddles

Adult paddle plates are removable to reveal pediatric paddles.

11130-000037



Internal Paddle Handles with Discharge Control

For use with LIFEPAK 15, 20e defibrillator/monitors.

11131-000001 (1 pair)



Internal handles/paddles

For use with LIFEPAK 20e. One-piece design includes both handle and paddle. 11131-000001 not compatible or required to use this item.

11131-000040 1.0 inch (1 pair)

11131-000041 1.6 inch (1 pair)

11131-000042 2.3 inch (1 pair)

11131-000043 3.0 inch (1 pair)

Therapy delivery accessories

EDGE System™ Electrodes for pacing/defibrillation/ECG with QUIK-COMBO® Connector

18 month minimum shelf life unless specified



EDGE System Electrodes with QUIK-COMBO Connector
11996-000091 (24" leadwire length)



EDGE System RTS (Radiotransparent) Electrodes with QUIK-COMBO Connector
11996-000090 (24" leadwire length)



Pediatric EDGE System RTS Electrodes with QUIK-COMBO Connector

For use only with manual monitor/defibrillators; 12 month minimum shelf life at time of shipment.
11996-000093 (24" leadwire length)



EDGE System Electrodes with QUIK-COMBO Connector and REDI-PAK® Preconnect System
11996-000017 (42" leadwire length)



QUIK-COMBO Therapy Cable
11110-000040

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

End-Tidal CO₂ (EtCO₂) monitoring accessories

Filterlines for intubated patients

Single patient use



FilterLine® SET

Key Applications: OR, EMS, ED, Rapid Response Teams, Transport.

Adult/Pediatric

11996-000081 (25/pk, 200 cm)
11996-000164 (25/pk, 400 cm) extra long
11996-000435 (100/pk, 200 cm)



FilterLine H SET

Key Applications: Critical care, Humidified Environments.

Adult/Pediatric

11996-000080 (25/pk, 200 cm)

Infant/Neonatal

11996-000001 (25/pk, 200 cm)

Non-intubated filterlines

Single patient use



Smart CapnoLine® Plus

Key applications: Procedural sedation, Lower GI procedures, MAC, EMS, ED, Rapid Response Teams.

Adult with O₂

11996-000163 (25/pk, 200 cm)
11996-000167 (100/pk, 200 cm)
11996-000165 (25/pk, 400 cm) extra long

Adult without O₂

11996-000162 (25/pk, 200 cm)
11996-000166 (100/pk, 200 cm)



Smart CapnoLine

Key applications: Procedural sedation, EMS, ED, Rapid Response Teams.

Pediatric with O₂

11996-000128 (25/pk, 200 cm)

Pediatric without O₂

11996-000120 (25/pk, 200 cm)

SpO₂ monitoring accessories


Masimo LNC® reusable patient cables



LNCS Patient Cable

11171-000024 (4 ft)
 11171-000016 (10 ft)
 11171-000025 (14 ft)
 11171-000027 (4 ft extension)

Masimo LNC compatible SpO₂ only sensors



LNCS® Adhesive Sensor

Single Patient Use. (20/box)
 11171-000019 (Ad)
 11171-000020 (Ped)
 11171-000029 (Neo/Pt)
 11171-000028 (Neo/Ad)
 11171-000031 (Inf)



LNCS® Reusable Sensor

11171-000017 (Ad)
 11171-000018 (Ped)



LNCS Reusable Soft Sensor

11171-000052 (Ad)

Masimo reusable direct connect sensors and cables




Direct Conect Reusable Soft Sensor

11171-000051 (Ad)

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

SpO₂ monitoring accessories

Additional accessories and adapter cables



**Reusable Ambient
Light Shield**
11171-000054 (5/pack)



**Disposable Ambient
Light Shield**
11171-000055 (10/pack)

SpO₂ monitoring accessories

Masimo to Nellcor™ adapter

To use any Nellcor accessories, your LIFEPAK 20e defibrillator/monitor must be configured for use with Nellcor technology at the time of LIFEPAK 20e defibrillator/monitor purchase.



MNC-1 Adapter Cable

Allows specially configured LIFEPAK 20e defibrillator/monitor with Masimo SpO₂ to connect to Nellcor sensors.

11996-000183 (10 ft)

11996-000198 (4 ft)

Masimo to Nellcor adapter compatible SpO₂ sensors



DURASENSOR Reusable Clip

11996-000060 (Adult - Ref# DS100A)

DURA-Y Multisite Reusable Sensor

11996-000106 (≥ 1 kg - Ref# D-YS)



Oxiband Reusable Sensor

Includes 50 disposable adhesive sensors.

11996-000061 (Ad/Neo)

11996-000062 (Ped/Inf)

Disposable Adhesive Bandage Wrap (100/pk)

Not Available in Canada.

11996-000048 (Ad/Neo - Ref# ADH A/N)

11996-000049 (Ped/Inf - Ref# ADH P/I)

Disposable Adhesive Foam Wrap

Not Available in Canada.

11996-000110 (Ad/Neo - Ref# FOAM A/N)

11996-000108 (Ped/Inf - Ref# FOAM P/I)



Oxisensor II (24/box)

Adult

11996-000113 (18 in. - Ref# D25)

11996-000114 (36 in. - Ref# D25L)

Pediatric

11996-000116 (18 in. - Ref# D20)



Oxisensor II (24/box)

Infant

11996-000115 (18 in. - Ref# I20)

Neonatal/Adult

11996-000117 (18 in. - Ref# N25)

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

Carrying cases and mounting options



Basic Carrying Case/Accessories Organizer

Includes set of right and left pouches for additional accessories.

11260-000018



Carrying Case for LIFEPAK 20e Defibrillator with Module

Includes a set of right and left pouches for additional storage of accessories. For use with a LIFEPAK 20e device with CodeManagement Module.

11260-000045



Top Pouch

For additional storage: insert in place of standard paddles.

11260-000043

Shoulder Strap for Basic Carrying Case

Includes mounting brackets and hardware.

11260-000041



Docking Station

Attaches to crash cart to lock LIFEPAK 20e defibrillator/monitor safely in place.

21330-000996

Right Side Accessory Pouch

Attaches to LIFEPAK 20e defibrillator/monitor and holds items such as QUIK-COMBO therapy cable, electrodes and SIGNAGEL electrode gel.

11260-000016

Training tools and testers




Patient Simulator (QUIK-COMBO)

Connects directly to your LIFEPAK defibrillator for safe simulation of cardioversion and electrical capture. Generates fibrillation, tachycardias and bradycardias, as well as ST segment and T wave abnormalities.

11996-000311 (12-Lead)

11996-000310 (3-Lead)



Defibrillator Checker

Tests integrity of energy delivery through the standard hard paddles. Neon light indicates energy has been delivered.

11998-000060



QUIK-COMBO Test Plug

Test plug connects to QUIK-COMBO therapy cable to test cable integrity.

11113-000002

LIFEPAK 20e defibrillator/monitor
with CodeManagement Module

Literature, communication and third party accessories



Inservice DVD
26500-003545 (English - DVD)

Service Manual
26500-003543 (English - CD - LIFEPAK 20e
with CodeManagement Module)

Operating Instructions:
LIFEPAK 20e Defibrillator/Monitor
with CodeManagement Module
26500-002570

Communication accessories



Transport Configuration Cable
For connecting LIFEPAK 20e to
LIFEPAK 20e defibrillator/monitor.
11230-000019 (6 ft)



LIFEPAK 20e



1
2
3
RES MODE
ANALYZE
LEAD SIZE SYNC
PACER
RATE
CURRENT
ALARMS
OPTIONS
PAUSE

PRINT CODE SUMMARY EVENT
AC Mains Service Speed Dial

ECG SpO2

WARNING: Possibility of electrical shock. For use only by qualified personnel. DANGER: Excessive heat. Do not use in the presence of flammable gases.

9

WARNING

**For further information, please contact Stryker at 800 442 1142 (U.S.),
800 668 8323 (Canada) or visit our website at strykeremergencycare.com**

Emergency Care

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CodeManagement Module, EDGE System, LIFEPAK, REDI-PAK, QUIK-COMBO, Stryker. Masimo, the Radical logo, Rainbow and SET are registered trademarks of Masimo Corporation. All other trademarks are trademarks of their respective owners or holders.

The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

GDR 3300292_K
Copyright © 2021 Stryker

Distributed by:

Stryker
11811 Willows Road NE
Redmond, WA, 98052 U.S.A.
Toll free 800 442 1142
strykeremergencycare.com

Distributed in Canada by:

Stryker Canada
2 Medicorum Place
Waterdown, Ontario
L8B 1W2
Canada
Toll free 800 668 8323