

stryker



Signature Portfolio

Maestro Air Motor Technology

The pneumatic power and feel you need

With the customization you want

As part of the Signature Portfolio, our latest Maestro Air motor reflects multiple refinements aimed at increasing its reliability and your overall experience. Rated up to 150 psi, today's motor provides the speed, power and torque to handle the most taxing drill scenarios, coupled with the intuitive feel and feedback you like from a pneumatic system. With a wealth of choices to suit your specialty and preferences, you can create a drill portfolio reflective of your own signature style.

Power and reliability

- Pneumatic motor rated up to 150 psi
- 20% more torque than previous model
- 25% more power in watts
- Speeds up to 75K RPMs

Customization

- Motors with or without handswitch
- Four foot pedal options
- The industry's most comprehensive array of attachments
- Industry-unique cutting accessories

Ergonomic freedom of motion

- Refined attachment design for improved balance and maneuverability
- Improved swivel at hose, increasing handpiece freedom of movement
- Lighter, more flexible hose with low-friction coating to reduce catching on surgical fabric and drapes, or collection of dust and debris



Dr. Lin
Custom preferences¹

Maestro Air
Universal foot pedal
130 psi
Irrigation rate 25%

Foot pedal options

Exhausting

Instantly vents pressurized gas at point of foot pedal, when foot is removed



Non-exhausting

Vents pressurized gas through motor system



Universal

Functions three ways in exhausting, non-exhausting and handswitch modes



Wide/heavy

For those who prefer wide-format foot controls



Compatible with the most expansive array of attachments and cutting accessories including:

- Multi-notch telescoping accessories (up to 10mm); reduces inventory up to 5:1
- 2Flute head design

Color coding instantly identifies matching attachments and cutting accessories

Alignment dots and Snap-Connect for easy and consistent assembly/disassembly

Design your ideal drill portfolio today

To learn more about Maestro Air motors, or the entire Signature Portfolio and its broad array of compatible attachments and cutting accessories, visit neurosurgical.stryker.com or call **800 253 3210**.

Ordering information

Part number	Product description
5407-200-000	Maestro Air Motor
5407-200-00H	Maestro Air Motor with handswitch
5400-207-000	Non-exhausting foot pedal
5400-208-000	Exhausting foot pedal
5400-218-000	Universal foot pedal
5400-207-000s1	Wide/heavy foot pedal
5400-200-002	DISS to Schrader air adapter
5400-200-003	DISS to ARO air adapter
5400-005-000	Maestro oil cartridge (10/pkg)
5400-006-000	Maestro diffuser (10/pkg)
5400-005-006	Oil cartridge/diffuser kit (10/pkg)
5400-206-000	Auto irrigator
5400-206-001	Irrigation cable

1. Surgeon names and preferences are for demonstration purposes only to help depict the range of performance and customizations possible.

Neurosurgical

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker Corporation or its affiliates own, use, or have applied for the following trademarks or service marks: Maestro and Stryker. All other trademarks are trademarks of their respective owners or holders.