

Pello

Take customization a step further

Designed with ENTs in mind

With one-of-a-kind-blades, Pello is designed for efficient cutting and one-click blade insertion and removal.¹

Balanced

The lightweight ergonomic design and comfortable pencil grip allows for precise positioning.¹

Reduced clogging

Pello's proprietary blade technology and internal shelf are designed to minimize clogging¹ for an improved surgical experience.



Pello

Blades and burs

Pello offers a variety of attachments, blades and burs designed to meet the needs of your everyday procedures.

Straight shavers



Pello 2.5mm straight shaver blade



Pello 3mm straight shaver blade



Pello 4mm straight shaver blade

Angled shavers



Pello 4mm angled 12 degree shaver blade




Pello 4mm angled 40 degree shaver blade



Pello 4mm angled 60 degree shaver blade

Burs



Pello 4mm angled 70 degree diamond bur



Pello 4mm angled 15 degree diamond bur

Pello

Powered by CORE 2

Stryker's CORE 2 console (Consolidated Operating Room Equipment) provides a convenient, single control center for a multitude of powered instruments and devices across multiple specialties with the ability to recall customized settings.²

The CORE 2 offers intuitive control with the following features:

Adjustable cutting speed

Easily increase or decrease RPM to match procedural requirements.

Selectable irrigation modes

Choose irrigation setting that best suit your needs.

Blade motion control

Switch between continuous and oscillating modes.

Custom setting recall

Save, store, and transfer user-defined presets for efficiency and consistency.



1. D0000019073,
2. CORE 2 console IFU 5400D0000220848-052-700, Pello IFU 700000890732

ENT

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery. The information presented is intended to demonstrate Stryker's products.

A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area. Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: CORE, Pello and Stryker. All other trademarks are trademarks of their respective owners or holders. The absence of a product, feature, or service name, or logo from this list does not constitute a waiver of Stryker's trademark or other intellectual property rights concerning that name or logo.