

Why **Hoffmann® LRF smartHEX™** matters




Patient compliance

Semi-automatic controllers combined with the mobile app notifications designed to support patient compliance. As the mobile app guides the controllers through the correction plan, the surgeon can see the progress made via the web app and support remote residual corrections when needed.

Fixation and instrumentation

Options of fixation components and instrumentation from tabbed rings, five sizes of double telescoping struts, HA coated pins, offset plates, wire posts and pin cubes, multiple pin adapters, clickers, wire guides, and more allow the surgeon to adapt to patient anatomy and correct severe deformities. Double telescoping struts increase length range and minimize the need for strut changeouts.



Tray and sterilization efficiency

Trays designed for Aesculap containers to organize and protect during sterilization. Added ability to pre-build frame prior to sterilization in accordance with the IFU.

Enhanced software

The smartHEX web application offers assistance in planning corrections for a range of deformities and complex limb reconstruction procedures. Line tool, complete joint angle analysis, overcorrection, and ring-to-ring distance computation designed to enable surgeons with flexibility in planning frame constructs and support deformity corrections.



Trauma & Extremities

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: HexConnect, Hoffmann, smartHEX, Stryker. All other trademarks are trademarks of their respective owners or holders.

This document is applicable to the United States of America.

FA-HLHX-OTHD-3861002_REV-0_en_us

Copyright © 2026 Stryker