

Hoffmann® LRF smartheX™



Bringing the future into the frame

Hoffmann® LRF smartheX™




Semi-automatic correction ability via Bluetooth enabled controllers that adjust the struts according to the surgeon's correction plan.

Mobile application downloads the correction plan for the patient and guides the controllers through the semi-automatic adjustments.




Stryker HexConnect



Additional components introduced to the system combine performance with simplicity: Engineered with high-performance materials machined to precise tolerances, LRF components are designed to be low profile and to offer enhanced physical access and radiographic visibility without compromising stability.

Reimagined software brings an all-new interface and powerful planning tools – from complete joint angle analysis to advanced complex deformity planning – so you can plan with confidence and efficiency.

To access the software, visit.stryker.com/smarthex-login
To register for an account, visit.stryker.com/insight-center



Tailored solutions offer the right construct for each case, with multiple configuration options and both manual and semi-autonomous correction plan capabilities.

Trauma & Extremities

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: HexConnect, Hoffmann, smartheX, Stryker. All other trademarks are trademarks of their respective owners or holders.

This document is applicable to the United States of America.

FA-HLHX-SS-3360650_REV-1_en_us
Copyright © 2026 Stryker