

Synchfix™ EVT

Syndesmotic fixation device

Plating compatibility guide

Synchfix EVT may be used standalone or with titanium plates designed to accept 3.5mm non-locking screws



Part #	Description
86SYN205	Synchfix EVT inserter with implant construct

Internal options

Ortholoc 2 Ankle Fracture LP


Syndesmosis holes are located in Ortholoc 3Di lateral fibula plates and the 4 hole syndesmosis plate

Syndesmosis and Lateral Fibula Plates

Part #	Description
5890KITB	Ortholoc 3Di Ankle FX Fibula Implant Kit
5890KITF	Ortholoc 3Di Ankle FX OFF Fibula Implant Kit

1. Straight syndesmosis plate
2. Offset syndesmosis plate
3. Offset lateral fibula plate
4. Straight lateral fibula plate

1.3mm thick



Pangea


Distal Lateral Fibula and Distal Posterolateral Fibula Plates

Part #	Description
542223	Distal Fibula Plate Tray

1. Distal lateral fibula plate
2. Distal posterolateral fibula plate



Distal cluster: 2mm thick
Shaft: 2.6mm thick




VariAx 2

Fibula and One Third Tubular Plates

Part #	Description
RSAM	VariAx 2 Small Fragment add-on tray

1. Distal lateral fibula plate
2. Straight fibula plate
3. One third tubular plate

Distal lateral fibula plate
Distal end: 1.3mm thick
Shaft: 2.0mm thick
Straight fibula plate
2.0mm thick
One third tubular plates
1.7mm thick



Elevating care with our next generation device
Experience. Vision. Technique.

Foot & Ankle

Visit us at footankle.stryker.com

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate a Stryker product. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any Stryker product. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use, or have applied for the following trademarks or service marks: Ortholoc, Pangea, Stryker, Synchfix, Synchfix EVT, VariAx. All other trademarks are trademarks of their respective owners or holders.

FA-GSNPS-SS-2663029_REV-0_en_us

Copyright © 2025 Stryker