



Novae cup instrument tray #1



No.	Catalog #	Description
	RA99990098	Novae Cup Instruments Tray 1
1	RA90370241 to RA90370269	Acetabular Reamer 41mm - 69mm
2	RA90670085	Reamer Handle
3	RA90270002	Reamer Adaptor
4	RA90430300	Clicker Liner Head Impactor
5	RA90560072	ICJB E Handle


No.	Catalog #	Description
6	RA90460171	Straight Gripper Shaft
7	RA90310182	Trial Cup Extractor
8	RA90250223	Cup Impactor
9	RA90250222	Liner Impactor
10	RAMP010324	CLICKER SOCLE

Novae cup instrument tray #2



No.	Catalog #	Description
1	RA90194322 to RA90196528	Retentive Trial Liner 43mm-65mm
2	RA90390943 to RA90390965	Trial Cup 43mm-65mm
3	RA90192228	Retentive Trial Liner Adaptor 22.2mm / 28mm
4	RA90460129	Minimally Invasive Gripper Shaft
5	RAMA030090	Guide Rod
6	RA90590143 to RA90590165	Cup Holder Adaptor 43mm-65mm
7	RA90310022	Liner Remover
8	RA90510007	Peg Forceps

Novae cup pegs and screw instruments tray



No.	Catalog #	Description
	RA99990062	Novae Pegs & Screw Instruments Tray
	RA99990106	Novae Pegs & Screw Instruments Lid
1	RA90310180	Peg Extractor / 6mm Hex Wrench / Flange Modeller
2	RA90430004	Straight Peg Impactor
3	RA90430020	Curved Peg Impactor
4	RA90450030	Depth Gauge

No.	Catalog #	Description
5	RA90710040	3.5mm Manual Hexagonal Screwdriver
6	RA90410037	Drill Guide for 3.2mm Dia.
7	RA90710035	3.5mm Motorised Hexagonal Screw Driver
8	RA90350001	3.2mm Dia. Screw Drill
9	RA90470022	5mm Dia. Flexible Drill Bit

A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker's global product offerings. A surgeon must always refer to the package insert, product label and/or instructions for use before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your sales representative if you have questions about the availability of products in your area.

Stryker Corporation or its affiliates own, use or have applied for the following trademarks or service marks: Novae, Stryker. All other trademarks are trademarks of their respective owners or holders.