


Triathlon[®]

Medial Stabilized Insert

A **trusted legacy** strengthened through **innovation**

Triathlon has a demonstrated history of clinical success¹⁻³, and we continue to build upon this legacy with innovative additions to the portfolio.

The **Triathlon Medial Stabilized Insert** is the latest of these additions, designed to deliver enhanced stability⁴⁻⁶ and kinematics⁷ to your Triathlon Total Knee procedures.



Experience the clinical advantages of our market-leading⁸ **Mako SmartRobotics™** and **Triathlon Cementless** technologies, now with your preferred tibial insert.

Triathlon Medial Stabilized Insert



Enhanced stability and kinematics

Triathlon's single radius design provides consistent collateral ligament tension and stability throughout active flexion.⁹⁻¹⁰ Designed with increased medial conformity,⁴⁻⁶ the Triathlon Medial Stabilized Insert builds upon this by delivering **enhanced stability**⁴⁻⁶ and more native rotational **kinematics**⁷ to your Triathlon primary knee patients, both with and without an intact PCL*.



Reproducible execution

Implant design and positioning work together to impact knee stability and kinematics.¹¹ Pair the Triathlon Medial Stabilized insert with patient specific **CT-based planning** and the precision¹² of **AccuStop™ haptic technology** of Mako SmartRobotics™. This combination allows for a personalized workflow, and helps you achieve your preferred alignment philosophy with confidence.



Trusted performance

The Triathlon MS Insert seamlessly joins the Triathlon Total Knee System, which has a proven history of delivering stability,¹⁰ satisfaction¹³ and excellent long-term survivorship.¹⁻³ This is highlighted by the **trusted fixation**^{14,15} of Triathlon Cementless. Multiple studies reported greater than 99% survivorship for aseptic loosening at 10 years^{14,15} and high patient satisfaction rates^{14,16} on top of its improved efficiency¹⁷⁻¹⁹ vs. cemented TKA.

*As demonstrated in controlled laboratory, computational and preclinical evaluations. Clinical outcomes have not been established.

References:

1. American Joint Replacement Registry (AJRR). 2025 annual report. Rosemont, IL: AAOS; 2025.
2. National Joint Registry. 22nd annual report 2025: Surgical data to 31 December 2024. 2025.
3. Australian Orthopaedic Association National Joint Replacement Registry (AOANJRR). Hip, knee & shoulder arthroplasty annual report. 2025.
4. Stryker internal report D0000373518.
5. Stryker internal report D0000361177.
6. Stryker internal report D0000373525.
7. Stryker internal report D0000367859.
8. Orthopedic Network News. Hip and knee implant review. 2025;36(3).
9. Cook L, et al. J Knee Surg. 2012;25:249-254.
10. Larsen L, et al. J Arthroplasty. 2015;30(6):1062-1067.
11. Angerame MR, et al. J Arthroplasty. 2019;34(10):2502-2510.
12. Deckey DG, et al. Bone Joint J. 2021;103-B(6 Suppl A):74-80.
13. Scott CEH, et al. Bone Jt Open. 2023;4(10):808-816.
14. Stoltz MJ, et al. J Arthroplasty. 2025; In press.
15. Zappley NR, et al. J Arthroplasty. 2025;40(7 Suppl 1):S250-S255.
16. Hannon CP, et al. J Bone Joint Surg Am. 2023;105(18):1430-1434.
17. Nam D, et al. J Bone Joint Surg Am. 2019;101(13):1185-1192.
18. Sharpe K, et al. Presented at ISTA; 2018; London, UK.
19. Lawrie CM, et al. Bone Joint J. 2019;101-B(7 Suppl C):61-63.

Joint Replacement

This document is intended solely for the use of healthcare professionals. A surgeon must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. We do not dispense medical advice and recommend that surgeons be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate Stryker's products. A surgeon must always refer to the package insert, product label and/or instructions for use, including the instructions for cleaning and sterilization (if applicable), before using any of Stryker's products. Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your representative if you have questions about the availability of Stryker's products in your area.

Stryker or its affiliated entities own, use, or have applied for the following trademarks or service marks: AccuStop, Mako, SmartRobotics, Stryker, Triathlon. All other trademarks are trademarks of their respective owners or holders.

JR-GSNPS-SS-2890650 Copyright © 2026 Stryker