

Furniture

Bedside Cabinet | Overbed Table

Complete your hospital room with the Furniture from Stryker

The Stryker Furniture consist of bedside cabinets, bedside cabinet with integrated overbed table and individual overbed tables. Different colour options allows the furniture to match with both the hospital bed and room environment.

Together with the SV2 hospital bed, Stryker can offer a complete and contemporary hospital room.

Bedside stands



MAC-307-S

Specifications

Weight	25 kg
Width	46 cm
Length	47 cm
Height	81 cm

Product features

- Modern stand with single drawer and cabinet
- Drawer with removable ABS plastic shelving tray
- Steel body with aluminum profile
- HPL top surface and front panel
- Cross-lockable castors for stability



MAC-308-S

Specifications

Weight	29 kg
Width	46 cm
Length	47 cm
Height	81 cm

Product features

- Modern stand with two drawers and cabinet
- Drawer with removable ABS plastic shelving tray
- Steel body with aluminum profile
- HPL top surface and front panel
- Cross-lockable castors for stability



MAC-310

Specifications

Weight	36.5 kg
Width	45 cm
Length	45.5 cm
Height	77 cm

Product features

- Integrated overbed table with two drawers
- Drawer with removable ABS plastic shelving tray
- Steel body with aluminum profile
- HPL top surface and front panel
- Cross-lockable castors for stability

Overbed tables



MOT-404

Specifications


Weight	15 kg
Width	42 cm
Length	80 cm
Height	80-104 cm

Product features

- Fixed top OBTable with U-Shaped base
- ABS plastic surface with raised edge for spill containment
- Gas spring mechanism for height adjustability
- Dual locable castors for stability
- Available in standard gray colour

Bedside Stands and Overbed Table

are available in 5 different colours:



MOT-406

Specifications

Weight	17 kg
Width	70 cm
Length	77 cm
Height	80-105 cm

Product features

- Single Side Tilt, Swivel Top OBTable with C-Shaped Base
- Gas spring mechanism for height adjustability
- Dual locable castors for stability
- Gray or white frame
- Available in five different colours

Medical

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.