



# Cub Paediatric Crib

Versatile. Access.  
Ease of use.



## When care matters most.

Responsible paediatric care demands a stretcher specifically designed around versatility, access and ease of use. The Cub Paediatric Crib is designed to support wellbeing-focused, efficient care practices while the warm, friendly appearance cultivates a pleasant environment for small children and their families.

### Swing-open access doors

allow quiet, convenient patient care without disturbing a sleeping child.

### Retracting backrest

helps eliminate head-end gap while raising the head end of the stretcher.

### Retractable fifth wheel

helps improve mobility.

### Counterbalance system

on siderail for easy raising and lowering.

### Activation points

All activation points are coloured red and require two adult size hand movements to initiate, helping ensure the patient remains in the crib.



### Retracting backrest

reduces the gap against the siderail while swing-open access doors allow easy patient care.



### Optional protective top

can be equipped with intermediate support position on each side.



### Lowered siderails

All four siderails can be lowered to provide unobstructed access which allows caregivers to have 360 degree access.

Multiple colour sticker options allow you to customise the stretcher for identification purposes.



# Cub Paediatric Crib

Versatile. Access. Ease of use.

## Standard features

- 10.2 cm (4") Ultra Comfort support surface
- Manually-operated pneumatic-assist backrest
- Manual foot prop
- Siderail with 22.9 cm (9") mandatory stop
- Counterbalance system on siderail for easy raising and lowering
- Oxygen bottle holder
- 15.2 cm (6") Omni Surface casters
- Central-locking brake system
- Dual side-mounted foot control hydraulics
- Uni-lower pedal
- Hands-free Trend./reverse Trend
- Retractable fifth-wheel steering
- Access doors at all four corners
- Fully-lowerable siderails provide unobstructed access around entire unit
- **Two year warranty**

## Optional features

- Siderail without 22.9 cm (9") mandatory stop
- Accessory supports with roller bumpers (head, foot or both ends)
- Chart shelf
- IV caddy
- Monitor tray
- Oxygen bottle retaining ring
- Siderail pads
- Retractable protective top with or without intermediate supports
- IV poles (permanent folding two-or three stage or standard removable)
- Upright oxygen bottle holder

## Specifications

<b>Model number</b>	FL19H
<b>Overall length &amp; width</b>	
No. accessory supports	166 x 94 cm (65.5" x 37")
Patient surface	184 x 95 cm (72.25" x 37.5")
<b>Weight capacity</b>	181.4 kg (400 lb.)
<b>Safe working load</b>	45.4 kg (100 lb.)
<b>Height range (to litter top)</b>	
High	101.6 cm (40")
Low	83.8 cm (33")
<b>Litter positioning</b>	
Manual activated	0° / 25° / 40° / 55°
Pneumatic assist	0° - 50°
Trend./reverse Trend.	±9.5°
Foot section angle	0° / 6° / 12°
<b>Support surface</b>	76 x 148 cm (30.13" x 58.25")
<b>Caster diameter</b>	15.2 cm (6")

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

For more information please contact your local sales representative.

**Stryker Australia Pty Ltd**  
8 Herbert Street  
St Leonards NSW 2065  
61 2 9467 1000  
www.stryker.com.au

**Stryker New Zealand Limited**  
PO Box 17-136  
Greenlane, Auckland, New Zealand  
64 9 573 1890  
www.stryker.com