

Emergency Relief Bed™

limited-release



Designed to meet the healthcare industry’s urgent and growing need for patient platforms during the current COVID-19 pandemic, our limited-release Emergency Relief Bed helps emergency responders and caregivers move and position patients efficiently.

Our Emergency Relief Bed is a low-cost, readily available solution intended to quickly serve those on the front lines of the COVID-19 crisis – ranging from triage centers and pop-up areas of care, to the emergency department, and all in-patient general hospital departments.

It features a 30-degree Head of Bed angle to accommodate patients in respiratory distress, including vented patients. It also features low height, a 6” thick mattress, and an attached IV pole.

Manufacturing target
10,000 beds per week

30° Head of Bed
Accommodates vented patients

6” Thick surface
Included

Part Number: 301900000000

Emergency Relief Bed Specifications

- Frame Dimensions: 86” x 36” x 11.5”
- 11.5” low height (without surface)
- 30° Head of Bed
- 5 Head of Bed angles
- Wheeled frame
- Caster lock brakes
- IV pole
- Weight capacity: 350 lbs.
- 6” thick foam mattress included

