

Fluoroscopy Stretcher



Access. Functionality. Convenience.

Fluoroscopy Stretcher's features are designed to simplify operation and reduce the need to reposition the patient during fluoroscopic imaging.



Pneumatic backrest

Fibers in pneumatic backrest at both ends improves imaging access and helps reduce the need to reposition a patient.



Enhanced imaging capabilities

U-shaped base and full-length C-arm access enhance imaging capabilities. High height litter option adds 18 cm to low and high height to accommodate different size C-Arms.



Patient access

Siderails fold down and out of the way for improved patient access.

Optional accessory rail accommodates surgical attachments.

4 wheel brake system helps secure treatment platform.

Fluoroscopy Stretcher



Standard features

- 2-inch (5 cm) mattress
- 6-inch (15 cm) nonconductive casters
- Dual end-mounted foot control hydraulics
- Dual fiberesin pneumatic backrest
- Fold-down siderails
- Four corner roller bumpers
- Four IV receptacles
- Four-wheel, steel-ring brake system with dual end activators
- Radiolucent litter top
- Locking-caster steering with dual end activators
- U-shaped base
- 2-stage permanent IV pole
- Integrated accessory rail

Optional features

- Footboard/chartholder
- IV poles
 - Permanent folding two- or three-stage at foot end
 - Removable
- Restraint straps
- Upright oxygen bottle holder

Specifications

Model number	1080
Overall length	87" (220.9 cm)
Overall width	
Siderails down	30.5" (77.4 cm)
Weight capacity	500 lb. (225 kg)
Height range*	
High	32.75" (83.2 cm)
Low	23.5" (59.7 cm)
Litter positioning	
Dual backrest operation	0°– 90°
Trend./Reverse Trend.	±8°
Patient surface	26" x 84" (66 x 213.4 cm)
Siderails	13" x 55" (33 x 139.5 cm)
Fluoroscopy window	27.5" x 66" (69.9 x 167.6 cm)
Caster diameter	6" (15 cm)
C-arm access width	54" (137.2 cm)

*High height litter option adds 7 inches to low and high height.

Medical

This document is intended solely for the use of healthcare professionals.

A healthcare professional must always rely on his or her own professional clinical judgment when deciding whether to use a particular product when treating a particular patient. Stryker does not dispense medical advice and recommends that healthcare professionals be trained in the use of any particular product before using it in surgery.

The information presented is intended to demonstrate the breadth of Stryker product offerings. A healthcare professional must always refer to the package insert, product label and/or instructions for use before using any Stryker product.

Products may not be available in all markets because product availability is subject to the regulatory and/or medical practices in individual markets. Please contact your Stryker representative if you have questions about the availability of Stryker products in your area.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holders.

The products depicted are CE marked in accordance with applicable EU Regulations and Directives.

This material is not intended for distribution outside the EU and EFTA.



01262005EN Rev F.3
DLS 12/2017
2017_16209