

Fluoroscopy Stretcher



Access. Functionality. Safety.

Our Fluoroscopy Stretcher's features are designed to simplify operation and minimize the need to reposition the patient during fluoroscopic imaging.



Pneumatic backrest

Fibers in pneumatic backrest at both ends improves imaging access and helps reduce the need to reposition a patient.



Enhanced imaging capabilities

U-shaped base and full-length C-arm access enhance imaging capabilities. High height litter option adds 7 inches to low and high height to accommodate different size C-Arms.



Patient access

Siderails fold down and out of the way for improved patient access.

Optional accessory rail accommodates surgical attachments.

4 wheel brake system helps secure treatment platform.

Fluoroscopy Stretcher



Standard features

- 2-inch (5 cm) mattress
- 6-inch (15 cm) nonconductive casters
- Dual end-mounted foot control hydraulics
- Dual fiberesin pneumatic backrest
- Fold-down siderails
- Four corner roller bumpers
- Four IV receptacles
- Four-wheel, steel-ring brake system with dual end activators
- Radiolucent litter top
- Locking-caster steering with dual end activators
- U-shaped base
- 2-stage permanent IV pole
- Integrated accessory rail

Optional features

- Footboard/chartholder
- IV poles
 - Permanent folding two- or three-stage at foot end
 - Removable
- Restraint straps
- Upright oxygen bottle holder

Specifications

Model number 1080

Overall length 87" (220.9 cm)

Overall width

Siderails down 30.5" (77.4 cm)

Weight capacity 500 lb. (225 kg)

Height range*

High 32.75" (83.2 cm)

Low 23.5" (59.7 cm)

Litter positioning

Dual backrest operation 0°– 90°

Trend./Reverse Trend. ±8°

Patient surface 26" x 84" (66 x 213.4 cm)

Siderails 13" x 55" (33 x 139.5 cm)

Fluoroscopy window 27.5" x 66" (69.9 x 167.6 cm)

Caster diameter 6" (15 cm)

C-arm access width 54" (137.2 cm)

*High height litter option adds 7 inches to low and high height.

Stryker reserves the right to change specifications without notice.

Stryker Corporation or its divisions or other corporate affiliated entities own, use or have applied for the following trademarks or service marks: Stryker. All other trademarks are trademarks of their respective owners or holder.

3800 E. Centre Avenue
Portage, MI 49002 USA
t: 269 329 2100
toll free: 800 327 0770

patienthandling.stryker.com

Copyright © 2017 Stryker
Mkt Lit-65 01262005 Rev E3
Printed in U.S.A.